

## 兔抗 ARF5 多克隆抗体

中文名称：兔抗 ARF5 多克隆抗体

英文名称：Anti-ARF5 rabbit polyclonal antibody

别名：ADP ribosylation factor 5

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ARF5

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene is a member of the human ADP-ribosylation factor (ARF) gene family. These genes encode small guanine nucleotide-binding proteins that stimulate the ADP-ribosyltransferase activity of cholera toxin and play a role in vesicular trafficking and as activators of phospholipase D. The gene products include 6 ARF proteins and 11 ARF-like proteins and constitute 1 family of the RAS superfamily. The ARF proteins are categorized as class I (ARF1, ARF2, and ARF3), class II (ARF4 and ARF5) and class III (ARF6). The members of each class share a common gene organization.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	ARF5
<b>Immunogen:</b>	Synthetic peptide of human ARF5
<b>Full name:</b>	ADP ribosylation factor 5
<b>SwissProt:</b>	P84085
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human colorectal cancer
<b>IHC Recommend dilution:</b>	30-150



