

兔抗 ARID1B 多克隆抗体

中文名称: 兔抗 ARID1B 多克隆抗体

英文名称: Anti-ARID1B rabbit polyclonal antibody

别名: OSA2; 6A3-5; DAN15; MRD12; P250R; BRIGHT; BAF250B; ELD/OSA1

储存: 冷冻 (-20℃)

抗原: ARID1B

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This locus encodes an AT-rich DNA interacting domain-containing protein. The encoded protein is a component of the SWI/SNF chromatin remodeling complex and may play a role in cell-cycle activation. The protein encoded by this locus is similar to AT-rich interactive domain-containing protein 1A. These two proteins function as alternative, mutually exclusive ARID-subunits of the SWI/SNF complex. The associated complexes play opposing roles. Alternatively spliced transcript variants encoding different isoforms have been described.
Applications:	ELISA, IHC
Name of antibody:	ARID1B
Immunogen:	Synthetic peptide of human ARID1B
Full name:	AT-rich interaction domain 1B

Synonyms:	OSA2; 6A3-5; DAN15; MRD12; P250R; BRIGHT; BAF250 B; ELD/OSA1
SwissProt:	Q8NFD5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100

