

## 兔抗 BBS12 多克隆抗体

中文名称: 兔抗 BBS12 多克隆抗体

英文名称: Anti-BBS12 rabbit polyclonal antibody

别名: Bardet-Biedl syndrome 12; C4orf24

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: BBS12

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is part of a complex that is involved in membrane trafficking. The encoded protein is a molecular chaperone that aids in protein folding upon ATP hydrolysis. This protein also plays a role in adipocyte differentiation. Defects in this gene are a cause of Bardet-Biedl syndrome type 12. Two transcript variants encoding the same protein have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BBS12
<b>Immunogen:</b>	Synthetic peptide of human BBS12
<b>Full name:</b>	Bardet-Biedl syndrome 12

<b>Synonyms:</b>	C4orf24
<b>SwissProt:</b>	Q6ZW61
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

