

## 兔抗 BBX 多克隆抗体

中文名称：兔抗 BBX 多克隆抗体

英文名称：Anti-BBX rabbit polyclonal antibody

别名：bobby sox homolog (Drosophila); HBP2; ARTC1; MDS001; HSPC339

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：BBX

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Transcription factor that is necessary for cell cycle progression from G1 to S phase.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	BBX
<b>Immunogen:</b>	Fusion protein of human BBX
<b>Full name:</b>	BBX high mobility group box domain containing
<b>Synonyms:</b>	HBP2; ARTC1; MDS001; HSPC339
<b>SwissProt:</b>	Q8WY36
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human prostate cancer
<b>IHC Recommend dilution:</b>	20-100

