

兔抗 CLDN2 多克隆抗体

中文名称：兔抗 CLDN2 多克隆抗体

英文名称：Anti-CLDN2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：CLDN2

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene product belongs to the claudin protein family whose members have been identified as major integral membrane proteins localized exclusively at tight junctions. Claudins are expressed in an organ-specific manner and regulate tissue-specific physiologic properties of tight junctions. This protein is expressed in the intestine. Alternatively spliced transcript variants with different 5' untranslated region have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	CLDN2
Immunogen:	Synthetic peptide of human CLDN2
Full name:	claudin 2

SwissProt:	P57739
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human gastric cancer
IHC Recommend dilution:	50-200



