

Article

Effects of Hormones and Epigenetic Regulation on the Callus and Adventitious Bud Induction of *Fraxinus mandshurica* Rupr.

Lei Yu ^{1,2}, Xinyu Li ^{1,2}, Hongmei Tian ³, Hualing Liu ⁴, Ying Xiao ⁴, Nansong Liang ^{1,2}, Xingtang Zhao ^{1,2} and Yaguang Zhan ^{1,2,*} 

¹ Key Laboratory of Saline-Alkaline Vegetation Ecology Restoration (SAVER), Ministry of Education, Northeast Forestry University, Harbin 150040, China; yuleibiology@126.com (L.Y.); lxy_3y6r@163.com (X.L.); liangnansong@nefu.edu.cn (N.L.); longhuxingsheng@163.com (X.Z.)

² Department of Forest Bioengineering, College of Life Science, Northeast Forestry University, Harbin 150040, China

³ Forest Botanical Garden of Heilongjiang Province, Harbin 150040, China; 18846830521@163.com

⁴ Shanhetun Forestry Bureau, Wuchang City, Harbin 150200, China; 18846831735@163.com (H.L.); xiaoyingsht@126.com (Y.X.)

* Correspondence: zhanyaguang2014@126.com; Tel.: +86-0451-82191753

Received: 26 April 2020; Accepted: 21 May 2020; Published: 24 May 2020



Abstract: *Fraxinus mandshurica* Rupr. (hereafter “*F. mandshurica*”) is known as one of northeast China’s important, valuable hardwood timber species. However, tissue culture and micropropagation of the species are difficult and have low efficiency, limiting asexual propagation. In this manuscript, stem explants were utilized to establish an effective regeneration system through adventitious bud organogenesis. The factors influencing callus regeneration *in vitro* were determined, and callus regeneration technology was established. The mechanism of adventitious bud formation was analyzed. Thidiazuron (TDZ) played a crucial role in the formation of adventitious buds. Elevated concentrations of TDZ were beneficial to callus induction and low concentrations of 6-benzyladenine (BA) led to loose state callus formation. The order of callus induction rates for different explants was stem cotyledon (100%) > segment (98.54%) > hypocotyl (92.56%) > root (50.71%). The effects of exogenous addition of 6-BA and TDZ on the endogenous hormone content of plants during the regeneration of adventitious buds were also assessed, as well as the expression characteristics of genes related to the regeneration pathway. The comprehensive analysis results showed that the suitable medium for callus induction and adventitious bud differentiation was c12 medium (MSB5 + 30 g/L sucrose + 7 g/L Agar + 5 mg/L 6-BA + 8 mg/L TDZ + 2 mg/L glycine + 0.1 mg/L IBA + 5% coconut water). The induction rates of callus and adventitious buds were 99.15% and 33.33%. The addition of 2.4 mg/L of the DNA demethylation reagent 5-azacytidine (5-aza) and 0.15 mg/L of the histone deacetylase inhibitor trichostatin A (TSA) increased the rates of adventitious bud induction by 17.78% over the control. This further laid the foundation for large-scale cultivation of excellent varieties and genetic transformation techniques.

Keywords: callus; adventitious shoot; regeneration; endogenous hormone contents; plant growth regulators; epigenetic regulation

1. Introduction

F. mandshurica Rupr. is an important broad-leaved tree species in northeastern China that can be used for various types of industrial and furniture timber and has significant economic and ecological value [1,2]. However, because of excessive deforestation, the number of adult *F. mandshurica*

declines every year. *F. mandshurica* has been listed as a rare species in China, and so, reducing the number of *F. mandshurica* removed by China's forestry development and promoting the breeding of *F. mandshurica*, have become imminent tasks [3]. *F. mandshurica* mainly relies on seed to propagate, but seeds have the characteristics of strong periodicity, deep dormancy, low germination rate, and slow propagation rate [2]. Seed germination must undergo a vernalization treatment for up to 6 months [4]. These factors increase the difficulty of *F. mandshurica* reproduction and limit its reproduction in forest development [5–7]. Asexual reproduction can maintain the excellent strains and genetic characteristics of *F. mandshurica*. Grafting, cutting and plant tissue culture can promote and apply fine varieties on a large scale. During the grafting process, the degree of joining between scions and the rootstock, the quality of the scions, and the method of grafting all affect the final survival rate of the plant. Because of these factors, it is difficult to gain ideal economic benefits [8]. Similarly, cutting is also subject to the above restrictions [9,10]. The propagation coefficient of *F. mandshurica* is often difficult to break through, so grafting and cutting cannot be used as the main means of asexual propagation. Genetic engineering technology provides an attractive way to effectively improve current varieties [11].

In plant tissue culture, exogenous hormones are applied to reprogram the cells to obtain regeneration ability. Firstly, callus formation is induced on call inducing medium (CIM) with a high auxin ratio, and then, on shoot inducing medium (SIM) with a high cytokinin ratio. *PLT2* and *PLT3* are two auxin signal response factors [12] that induce callus formation in Arabidopsis. Wound treatment can promote callus formation in all parts of Arabidopsis seedlings by inducing the expression of the AP2-ERF transcription factor *WIND1* [13,14]. Wound treatment can activate the expression of *YUC1* and increase the level of auxin [15]. Auxin can activate the expression of the WUSCHEL-related gene *WOX*. This mechanism can maintain a high level of auxin in cells with regeneration ability. In the process of bud regeneration, the signal responses of endogenous auxin and cytokinin are antagonistic to each other. The transcription factor *WUS* plays a key role in Arabidopsis bud regeneration [16]. The cytokinin pathways *ARR2* and *REV* up-regulate the expression of the *WUS* gene [17].

In epigenetic studies, histone acetylation is associated with auxin-mediated somatic embryogenesis in Arabidopsis [18,19]. Low DNA methylation levels favor Arabidopsis regeneration, and histone acetylation is associated with auxin-mediated Arabidopsis somatic cell developmental embryo development programs [18,20]. This *F. mandshurica* functional gene research can only investigate gene cloning and expression analysis, lacking effective authentication on the function of the method [21,22]. Concerning asexual propagation of *F. mandshurica*, the differentiation of adventitious buds in the callus pathway and the complete establishment of a regeneration system have not yet been reported. Therefore, establishing efficient and effective callus induction and an adventitious bud regeneration system in *F. mandshurica* is very important. Therefore, it is of profound significance to analyze the breeding and gene function of *F. mandshurica* through callus regeneration.

2. Materials and Methods

2.1. Disinfection and Inoculation of *F. mandshurica* Seeds

The explant material was chosen from the hypocotyls, roots and cotyledons of the sterile seedlings of *F. mandshurica* and the stem sections of the sterile tissue culture seedlings of the bristles cultured in this laboratory. Fully-fledged seed of *F. mandshurica* was selected and treated with winged sand at 4 °C for 3 months. The treated seeds were deprived of the outer seed coat and rinsed with running water. They were immersed in 70% alcohol for 30–45 s, and 10% NaClO for 15–20 min. They were rinsed four times on a sterile table with sterile water, the embryo was peeled off, and 20 g/L sucrose and 5.3 g/L agar were added. They were placed in hormone-containing WPM induction medium, and the infection rate and initiation rate of different disinfection combinations were recorded. Ten tissue materials were treated each time, and the experiments were repeated three times. Data were analyzed using SPSS 23.0 statistical software (IBM-SPSS 2015). The mean with standard error (\pm SE) is presented. The percent of callus formation, shoot bud induction, and number of shoots were subjected to analysis

of variance (ANOVA). Significant difference between treatments was tested by Duncan's multiple comparison test ($p = 0.05$).

2.2. Adventitious Bud Induction of Hypocotyl, Root, Cotyledon and Stem Segments of *F. mandshurica*

The aseptic seedlings that had been germinated for 21 days were selected, and the roots of the embryos were cut off. The hypocotyls, roots and cotyledons were inoculated into separate hormone combination media. Tissue culture that was growing well was selected, and the stem segment was inoculated into the culture medium between the two buds. The specific composition is shown in Table 1. The adventitious buds were formed by light after dark culture for two weeks. The cultures were incubated in growth chamber with 16:8 h light/dark period. The temperature was adjusted at 18 ± 2 °C for dark and 25 ± 2 °C for light periods. Ten tissue materials were treated each time, and the experiments were repeated three times. The callus induction rate at 20 days and the adventitious bud induction rate at 30 days were counted. In the induction process, the 5-aza to c12 medium was added to inhibit the DNA methylation level of the material and observe its effect on the adventitious bud induction rate. Auxin and cell division after the formation of adventitious buds in the stem segments of *F. mandshurica* culture tissue culture without adding 5-aza (c12 medium) and adding 2.4 mg/L 5-aza for 30 d were also further assessed.

Table 1. Adventitious bud induction medium of stem segments.

Media Number	Basic Medium Type	Sucrose (g/L)	Agar (g/L)	6-BA (mg/L)	TDZ (mg/L)	IBA (mg/L)	Coconut (v/v)	TSA (mg/L)	5-aza (mg/L)
c1	MSB5	30	7	0	6	0.1	5%	0	0
c2	MSB5	30	7	5	6	0.1	5%	0	0
c3	MSB5	30	7	6	6	0.1	5%	0	0
c4	MSB5	30	7	7	6	0.1	5%	0	0
c5	MSB5	30	7	8	6	0.1	5%	0	0
c6	MSB5	30	7	0	7	0.1	5%	0	0
c7	MSB5	30	7	5	7	0.1	5%	0	0
c8	MSB5	30	7	6	7	0.1	5%	0	0
c9	MSB5	30	7	7	7	0.1	5%	0	0
c10	MSB5	30	7	8	7	0.1	5%	0	0
c11	MSB5	30	7	0	8	0.1	5%	0	0
c12	MSB5	30	7	5	8	0.1	5%	0	0
c13	MSB5	30	7	6	8	0.1	5%	0	0
c14	MSB5	30	7	7	8	0.1	5%	0	0
c15	MSB5	30	7	8	8	0.1	5%	0	0
b12	MS	30	7	5	8	0.1	5%	0	0
d12	WPM	30	7	5	8	0.1	5%	0	0
T5	MSB5	30	7	5	8	0.1	5%	0.15	0
T6	MSB5	30	7	5	8	0.1	5%	0.3	0
T7	MSB5	30	7	5	8	0.1	5%	0.6	0
T8	MSB5	30	7	5	8	0.1	5%	1.2	0
A5	MSB5	30	7	5	8	0.1	5%	0	2.4
A6	MSB5	30	7	5	8	0.1	5%	0	6.1
A7	MSB5	30	7	5	8	0.1	5%	0	12.21
A8	MSB5	30	7	5	8	0.1	5%	0	24.42

The callus and adventitious bud medium was composed of MSB5 + 30 g/L sucrose + 7 g/L agar + 5 mg/L 6-BA + 8 mg/L TDZ + 2 mg/L-glycine + 0.1 mg/L IBA + 5% coconut water + 0.15 mg/L TSA. The secondary proliferation medium was composed of WPM + 20 g/L sucrose + 7 g/L agar + 8 mg/L 6-BA. The rooting medium was composed of WPM + 20 g/L sucrose + 7 g/L agar + 1.4 mg/L IBA + 0.7 mg/L NAA. After domestication, they were transplanted to a non-woven bag containing peat soil:vermiculite 3:1 mixed matrix. The seedling box was placed in culture conditions at a temperature of 25 °C, humidity 90%, light intensity $80 \mu\text{mol}\cdot\text{m}^{-2}\cdot\text{s}^{-1}$, light-dark alternate culture in a greenhouse with a light time of 16 h.

2.3. Extraction and Analysis of Endogenous Hormones in Plants

Samples were ground to a fine powder in liquid nitrogen with a mortar and pestle, and suspended in potassium phosphate buffer (50 mM, pH = 7.4). After vigorous shaking, they were centrifuged for 30 min at 8000 rpm at 4 °C. Endogenous hormones were extracted from plant extract samples using Milbio enzyme-linked bio-plant gibberellin (GA₃), zeatin nucleoside (ZR), abscisic acid (ABA) and indole acetic acid (IAA) enzyme-linked immunosorbent assay (ELISA) kits. For the determination, three biological replicates were set for each sampling point and three technical replicates were set. Specific methods were as follows: loading, enzyme, incubated with a solution, washed, color, termination, and measurement.

2.4. DNA Demethylation and Histone Acetylation

The total DNA was extracted from the leaves of *F. mandshurica* by CTAB method. The overall methylation level of DNA was determined using a ZYMO RESEARCH 5-mC DNA ELISA Kit (D5325). Three biological replicates were set for each sampling point and three technical replicates were set. The specific experimental methods were as follows: DNA coating, blocking, adding antibody and incubating at 37 °C for 1 h, color development for 30 min, and ELISA plate reading to measure the absorbance at 405–450 nm.

Histone deacetylase HDAC has the effects of regulating cell proliferation, inducing cell differentiation and apoptosis. In this experiment, the histone deacetylase inhibitor trichostatin A (TSA) was used to cultivate adventitious buds from callus stems of the *F. mandshurica*. HDAC enzyme activity gradually decreased, and histone acetylation level gradually increased. TSA can advance the induction of adventitious shoots. Enzyme activity assays were performed on the extracted proteins using the BioVision HDAC Active Fluorescence Assay Kit. The specific experimental methods were as follows: 1. Dilute the test sample (10–50 µg nuclear extract or cell lysate) to 85 µL (final volume) with ddH₂O in each well (only 85 µL ddH₂O is added for the blank control), and the positive control is diluted with 83 µL ddH₂O with 2 µL of HeLa nuclear extract, and the negative control had 83 µL of ddH₂O with 2 µL of trichostatin A or a sample containing no HDAC activity. 2. Add 10 µL of 10× HDAC Assay Buffer. 3. Add 5 µL of HDAC Fluorescent Substrate and mix well. 4. Incubate for 30 min at 37 °C. 5. Add 10 µL of lysine developer and incubate for 30 min at 37 °C. 6. Read the sample in a fluorescent plate reader (Ex/Em = 350–380 nm/440–460 nm). Histone deacetylase activity is expressed as the relative fluorescence value of each sample protein.

2.5. Quantitative Real-Time PCR Analyses

The RNA of the sample was extracted by the Tris-CTAB method, and the RNA was reverse transcribed into cDNA by TransScript All-in-One First-Strand cDNA Synthesis SuperMix for qPCR (One-Step gDNA Removal), and bioedited from the transcriptome data of *F. mandshurica*. The fluorescent quantitative PCR primers were designed according to the gene sequence (Table 2) and the quantitative amplification of the Applied Biosystems ABI 7500 quantitative PCR instrument was performed by TransStart Top Green qPCR SuperMix (+Dye II). The reaction system included template cDNA 1 µL (10 ng·µL⁻¹), 2× SYBR Premix Ex taq 10 µL, Rox reference DyeII 0.4 µL, primer F 0.8 µL, primer R 0.8 µL, and ddH₂O, added to 20 µL. The reaction procedure was: 95 °C for 1 min; 95 °C for 10 s, 60 °C for 30 s, 72 °C for 35 s, for 40 cycles. The tests were repeated three times per sample. The expression level was calculated by the 2^{-ΔΔCT} method using the *F. mandshurica* vascular protein gene *β-tubulin* as an internal reference gene. The experimental data were analyzed by SPSS23 software.

Table 2. Real-time PCR primers.

Primer Name	Primer Sequence
qFmTU-F	5'-AGGACGCTGCCAACAACCTTT-3'
qFmTU-R	5'-TTGAGGGGAAGGGTAAATAGTG-3'
qFmAHP6-F	5'-GTAAGTCACCACCAATGCCCA-3'
qFmAHP6-R	5'-TCAGAGTTACAGAGGGTGATGCC-3'
qFmYUC1-F	5'-CCAGCATCACCACCTTCAACAT-3'
qFmYUC1-R	5'-CGAACAGACCCTCGGCAAGTA-3'
qFmPLT2-F	5'-TGACACTCCGCTTCATAACCA-3'
qFmPLT2-R	5'-AACAGATAATCGGAACCACTCG-3'
qFmPLT3-F	5'-GGAACATTTGCTACTGAGGAGG-3'
qFmPLT3-R	5'-GAGATTATGATGATAGAATGGGGC-3'
qFmARR1-F	5'-GGGTTTGATGAGGAGAGTGGG-3'
qFmARR1-R	5'-GTCCCATTTAGTCCCCGTTT-3'
qFmARR2-F	5'-CTATGCCTGGGATGACTGGA-3'
qFmARR2-R	5'-GGAAAGCCATGATCGTGA-3'
qFmESR2-F	5'-GGCTTGGAACCTTTGACACG-3'
qFmESR2-R	5'-TGAAGCTATACGAACGGACAAG-3'
qFmPHV-F	5'-CGGATAAGCAAACTGGGGGT-3'
qFmPHV-R	5'-AGGGTGGTGGCTCCATTATTCA-3'
qFmREV-F	5'-CGGTTACTTTCCGCCAGCATA-3'
qFmREV-R	5'-GGCTTCATTCCAGGCATTTG-3'
qFmWIND1-F	5'-TGCCAATCACAGGATCTAAA-3'
qFmWIND1-R	5'-CTGACCAAAGTACACGAAA-3'
qFmWOX4-F	5'-TGGATTCTGGGAGCAACATA-3'
qFmWOX4-R	5'-AGCCGATTCCTTCTTGAT-3'
qFmWOX5-F	5'-GGCAATAATGGCGTCAGTGG-3'
qFmWOX5-R	5'-CCCTGAATGAGTTGAGTGGAAATA-3'
qFmWUS-F	5'-GGAGGAAGTCTTTGTGGTGC-3'
qFmWUS-R	5'-GAACACTTTGATTGGAGGTCC-3'
qFmZPR3-F	5'-GCAGTGCTTGATTCTCCTGGTT-3'
qFmZPR3-R	5'-ATGGAGAGGCTAAACTCAAAGCTG-3'

3. Results

3.1. Effect of Exogenous 6-BA and TDZ Hormone Combination on Adventitious Bud Induction

The effects of different combinations of 6-BA and TDZ hormones on callus and adventitious buds on the use of hypocotyls, cotyledons, roots and stems as explants were explored (Figure 1). The results of 20 hormone treatments on cotyledons and roots used as explants are shown in Table 3. Under a4 medium (Figure 1d), the cotyledon callus induction rate reached 100%, and the texture of the cotyledon callus was hard. The color was dark green and the outer layer was a light green. Under a20 medium (Figure 1c), the highest callus induction rate of the roots was 50.71%. The callus formed by the roots had a velvety texture with canary yellow color and severe browning of the internal callus. Adventitious buds were not induced in all combinations of the two explants.

The results of 22 hormone treatments on hypocotyl used as explants are shown in Table 4. The callus induction rates in the a12 (Figure 1a) and a20 (Figure 1b) media were 93.56% and 88.69%. Germination rates of the a4 (Figure 1e) and a20 (Figure 1f) media were 3.13% and 10.52%. When the concentration of 6-BA was the same, with the increase in TDZ concentration, the callus induction rate also increased, which shows that TDZ played an important role in the induction process of hypocotyl callus. The texture of callus formed by low-concentration 6-BA-induced material was densely opaque green, and the texture surface of callus formed by high-concentration 6-BA-induced material was relatively loose and emerald green. The higher the concentration, the looser the callus. Therefore, a20 medium is suitable as a hypocotyl callus induction and adventitious bud regeneration medium.

Table 3. Effect of hormone treatment on cotyledon and root callus induction and adventitious bud induction.

Media Number	6-BA (mg/L)	TDZ (mg/L)	Cotyledon Survival Rate (%)	Cotyledon Callus Induction Rate (%)	Root Survival Rate (%)	Root Callus Induction Rate (%)	Adventitious Bud Induction Rate (%)
a1	0	0	17.5 ± 5.46hi	55.00 ± 7.07d	25.42 ± 4.28jkl	0g	0
a2	0	0.1	50 ± 1.57f	0g	28.48 ± 5.07ijk	14.58 ± 7.92e	0
a3	0	0.5	59.52 ± 0.49ef	41.67 ± 3.27e	54.94 ± 1.42bcd	23.33 ± 2.35d	0
a4	0	1	100a	100a	57.71 ± 2.96abc	48.49 ± 1.44a	0
a5	1.5	0	12.5 ± 8.92i	0g	18.06 ± 0.39l	0g	0
a6	1.5	0.1	88.89 ± 8.17b	0g	47.17 ± 4.59def	10.42 ± 7.5ef	0
a7	1.5	0.5	55.56 ± 6.29f	44.44 ± 1.02e	51.39 ± 8.58cde	33.68 ± 4.15c	0
a8	1.5	1	78.42 ± 3.75cd	92.71 ± 3.15b	61.86 ± 0.69ab	27.42 ± 5.54cd	0
a9	2	0	27.27 ± 1.31gh	18.18 ± 6.96f	64.44 ± 6.81a	0g	0
a10	2	0.1	57.14 ± 9.24f	0g	36.31 ± 2.2gh	4.17 ± 7.48fg	0
a11	2	0.5	75.25 ± 6.75d	63.89 ± 8.46c	48.33 ± 7.63def	14.17 ± 2.68e	0
a12	2	1	87.5 ± 6.5bc	91.67 ± 6.29b	44.55 ± 6.8efg	40.97 ± 0.02b	0
a13	3	0	50 ± 8.05f	0g	28.65 ± 3.5ijk	0g	0
a14	3	0.1	20 ± 5.19hi	12.5 ± 3.33f	24.64 ± 0.83jkl	0g	0
a15	3	0.5	10 ± 6.44i	0g	33.93 ± 0.16hij	0g	0
a16	3	1	70.00 ± 3.57d	70 ± 5.11c	33.22 ± 0.06hij	29.64 ± 0.85cd	0
a17	5	0	32.14 ± 8.89g	0g	20.35 ± 6.81kl	0g	0
a18	5	0.1	33.18 ± 3.79g	15.48 ± 4.83f	40.24 ± 0.08fgh	3.33 ± 9.79fg	0
a19	5	0.5	68.75 ± 7.25de	48.75 ± 7.91de	27.56 ± 8.77ijk	10 ± 7.42ef	0
a20	5	1	50 ± 3.17f	16.67 ± 4.99f	58.47 ± 8.43abc	50.71 ± 0.05a	0

Mean ± SE followed by the same letter in the same column were not significantly different by Duncan's multiple comparison test ($p < 0.05$).



Figure 1. Hypocotyl, root and cotyledon-induced callus and hypocotyl regeneration buds. (a) Hypocotyls induced callus in a12 medium for 20 d (scale 5 mm); (b) Hypocotyls induced callus in a20 medium for 20 d; (c) Root callus induced in a20 medium; (d) Cotyledons induced callus in a4 medium; (e) Hypocotyls induced adventitious buds in a4 medium for 30 d (scale 5 mm); (f) Hypocotyls induced adventitious buds in a20 medium for 30 d (scale 2 mm).

Table 4. Effects of hormone treatment on hypocotyl callus and adventitious bud induction.

Media Number	6-BA (mg/L)	TDZ (mg/L)	Survival Rate (%)	Callus Induction Rate (%)	Adventitious Bud Induction Rate (%)
a1	0	0	26.39 ± 7.01de	10 ± 1.2ij	0c
a2	0	0.1	30.95 ± 3.58de	25 ± 8.41gh	0c
a3	0	0.5	82.59 ± 8.48b	77.38 ± 4.81b	0c
a4	0	1	96.43 ± 1.39a	71.43 ± 9.92b	3.13 ± 1.25b
a5	1.5	0	0f	0j	0c
a6	1.5	0.1	25 ± 6.83e	25 ± 6.36gh	0c
a7	1.5	0.5	58.33 ± 3.05c	53.33 ± 3.1cd	0c
a8	1.5	1	73.86 ± 7.51b	87.5 ± 6.24a	0c
a9	2	0	0f	0j	0c
a10	2	0.1	36.36 ± 6.73d	18.18 ± 4.42hi	0c
a11	2	0.5	78.47 ± 9.12b	72.69 ± 2.8b	0c
a12	2	1	95.83 ± 1.77a	93.56 ± 2.39a	0c
a13	3	0	28.57 ± 0.72de	37.5 ± 4.71ef	0c
a14	3	0.1	52.08 ± 9.97c	31.94 ± 1.81fg	0c
a15	3	0.5	30.16 ± 9.87de	0j	0c
a16	3	1	62.5 ± 2.03c	59.29 ± 0.84c	0c
a17	5	0	30 ± 7.14de	20 ± 7.48h	0c
a18	5	0.1	52.05 ± 2.08c	27.56 ± 7.24gh	0c
a19	5	0.5	78.61 ± 3.54b	61.67 ± 9.79c	0c
a20	5	1	81.75 ± 5.06b	88.69 ± 9.64a	10.52 ± 2.12a
a21	5	2	28.64 ± 6.22c	46.35 ± 1.67de	0c
a22	5	2	60.34 ± 7.93c	40.72 ± 8.65ef	0c

Mean ± SE followed by the same letter in the same column were not significantly different by Duncan's multiple comparison test ($p < 0.05$).

The results of 15 hormone treatments on stem sections of tissue culture seedlings as explants are shown in Table 5. The callus induced in the c12 medium was hard inside, and visible protuberances and multiple outgrowths, which subsequently developed into adventitious shoot buds, were observed (Figure 2a–d). The callus induction rate of stem explants in c12 medium was 99.15%, and the bud induction rate was 33.33%. Although callus induction was efficient, the adventitious bud differentiation rate was still not ideal.

Table 5. Effect of hormone treatment on stem callus and adventitious bud induction.

Media Number	6-BA (mg/L)	TDZ (mg/L)	Survival Rate (%)	Callus Induction Rate (%)	Adventitious Bud Induction Rate (%)
c1	0	6	35.75 ± 1.02f	67.98 ± 0.33k	0g
c2	5	6	65.74 ± 3.08cd	87.45 ± 0.22h	0g
c3	6	6	78.68 ± 6.79a	90.41 ± 0.07f	0g
c4	7	6	67.48 ± 0.44bcd	91.34 ± 0.45e	4.55 ± 0.15ef
c5	8	6	79.64 ± 4.3a	96.54 ± 0.7c	3.33 ± 1.23f
c6	0	7	45.67 ± 6.27ef	74.12 ± 0.84j	6.25 ± 2.14de
c7	5	7	78.14 ± 7.54ab	98.54 ± 0.41a	3.7 ± 1.33f
c8	6	7	77.64 ± 5.69ab	97.35 ± 0.09b	0g
c9	7	7	64.12 ± 4cd	88.97 ± 0.85g	11.81 ± 1.49c
c10	8	7	61.47 ± 1.59d	98.74 ± 0.37a	16.21 ± 1.19b
c11	0	8	41.24 ± 8.43ef	75.64 ± 0.03i	6.94 ± 1.5d
c12	5	8	74.65 ± 8.32abc	99.15 ± 0.27a	33.33 ± 1.17a
c13	6	8	74.51 ± 3.69abc	97.64 ± 0.36b	10.48 ± 1.45c
c14	7	8	66.66 ± 9.88cd	98.42 ± 0.03a	2.86 ± 1.6f
c15	8	8	50.79 ± 5.99e	95.11 ± 0.36d	9.72 ± 1.69c

Mean ± SE followed by the same letter in the same column were not significantly different by Duncan's multiple comparison test ($p < 0.05$).




Figure 2. Stem segment induced callus and its regenerated shoots. (a) The callus formed by the epidermis of the stem segment was cultured in C12 for 20 d (scale 5 mm); (b) Cultured for 30 d; (c) The callus formed by the wound of the stem segment was cultured in C12 for 20 d (scale 2 mm). (d) Regenerated shoots formed in 30 d (scale 5 mm).

3.2. Effects of 5-aza and TSA on Adventitious Bud Differentiation of *F. mandshurica*

DNA 5-mC content was detected and calculated by the DNA methylation detection kit. The results showed that the addition of 2.4 mg/L 5-aza was most beneficial to the induction of adventitious buds (Figure 3), and the degree of DNA methylation was the lowest (Figure 5a). The DNA methylation level was related to the 5-mC content, where a higher 5-mC content indicates a higher DNA methylation level (Figure 5b). The induction rate of adventitious buds was 50.79%, which was 17.46% higher than that without 5-aza. The callus induced by the stem segment was mostly yellow-green with no signs of aging, and the dense induction of adventitious buds had a significant effect. The expression of genes related to regeneration was analyzed. It was found that the cytokinin pathway gene *FmREV*, the shoot tip and root tip gene *FmWOX4*, and the auxin synthesis gene *FmYUC1* were significantly increased during regeneration after the addition of 5-aza, which were 1.30, 2.03, 3.33 times, respectively, that of the control (Figure 5c).

The results showed that TSA could inhibit HDAC enzyme activity (Figure 4) and promote the induction of adventitious buds (Figure 5d,e). The regeneration rate of adventitious buds induced by the addition of 0.15 mg/L TSA reached a maximum of 51.11%, which was 17.78% higher than without TSA. The callus induced was emerald green, and the adventitious buds were fairly young. The expression levels of *FmPLT3*, *FmREV*, *FmWOX5*, and *FmZPR3* genes were significantly increased after adding TSA, which were 1.64, 8.82, 1.4 and 1.32 times higher than those of the control group (Figure 5f). In conclusion, adding 2.4 mg/L 5-aza or 0.15 mg/L TSA contributed to the induction of adventitious bud differentiation and the dense adventitious buds in the callus of tissue culture seedlings.




Figure 3. Adding different concentrations of 5-aza to induce the formation of adventitious buds. (a) Adding 2.4 mg/L 5-aza; (b) Adding 6.1 mg/L 5-aza; (c) Adding 12.2 mg/L 5-aza; (d) Adding 24.42 mg/L 5-aza (scale 2 mm).




Figure 4. Adventitious buds induced by the addition of different concentrations of TSA. (a) Adding 0.15 mg/L TSA; (b) Adding 0.3 mg/L TSA; (c) Adding 0.6 mg/L TSA; (d) Adding 1.2 mg/L TSA (scale is 2 mm).




Figure 5. Epigenetic and related gene expression levels during stem bud regeneration. (a) The 5-aza concentration affects the 5-mC level of shoot bud regeneration during 30 d of culture. (b) Changes in DNA methylation levels at various stages of regeneration of the *Fraxinus mandshurica* Rupr. (c) 5-aza affects the expression level of regenerative bud-related genes in stem segments. (d) The TSA concentration affects the level of regenerated bud HDAC in shoots of stem segments during 30 d of culture. (e) TSA affects the change in HDAC activity during regeneration. (f) TSA affects the expression level of regenerative bud-related genes in stem segments. Different letters above bars within statistically significant differences between different times of the treatments at the $p < 0.05$ level according to Duncan's multiple range test.

3.3. Effects of TDZ and 6-BA on Endogenous Hormone Levels in the Regeneration of *F. mandshurica*

After adding TDZ externally, with the increase in TDZ concentration, the content of IAA, ZR and ABA increased (Figure 6a–d). The content of GA₃ first increased and then decreased. IAA/ABA, ZR/ABA and GA₃/ABA content ratios increased, and ZR/IAA content ratios decreased (Figure 6e–h). It can be seen that exogenous application of 8 mg TDZ can promote the synthesis or accumulation of IAA, Zr and GA₃, and inhibit the synthesis or accumulation of ABA. Compared with the control group, IAA accumulation increased 1.69 times from the original level, GA₃ increased 1.33 times, ZR increased 1.29 times and ABA decreased 0.48 times. Analysis, combined with Table 5, showed that the change of adventitious bud induction rate was consistent with the change of IAA, ZR, and GA₃ hormone levels, indicating that TDZ can increase the levels of IAA, ZR, GA₃ hormones and regulate the differentiation of adventitious buds.




Figure 6. Effect of TDZ on endogenous hormone levels in regeneration of *F. mandshurica*. (a–d) CK was the stem tissue, and C1, C6 and C11 were the three endogenous hormones measured under the same culture conditions, and the TDZ concentrations were 6 mg/L, 7 mg/L and 8 mg/L for 30 d. (e–h) The relative amount of endogenous hormones. Different letters above bars within statistically significant differences between different times of the treatments at the $p < 0.05$ level according to Duncan's multiple range test.

As the concentration of exogenous 6-BA gradually increases, the content of endogenous IAA decreases, the content of ZR decreases slightly, the content of ABA increases, and the content of GA₃ remains stable (Figure 7a–d). The ratios of IAA/ABA, ZR/ABA and GA₃/ABA also reduced. The ratio of ZR/IAA content showed a “W”-shaped fluctuation with the increase in 6-BA concentration (Figure 7e–h). A low concentration of 5–7 mg/L 6-BA stimulated the synthesis or accumulation of IAA, ZR and GA₃, and inhibited the synthesis or accumulation of ABA in the 6-BA range. According to the analysis of the adventitious bud induction rate in Table 5, it was found that the adventitious bud induction rate was consistent with the change in IAA hormone levels, indicating that exogenous 6-BA plays a role in regulating the endogenous hormones IAA and ZR.




Figure 7. Effect of 6-BA on endogenous hormone levels in regeneration of *F. mandshurica*. (a–d) CK was the stem tissue, and C11, C12, C13 and C14 were the four endogenous hormones measured under the same culture conditions, and the 6-BA concentrations were 0, 5, 6 and 7 mg/L for 30 d. (e–h) The relative amount of endogenous hormones. Different letters above bars within statistically significant differences between different times of the treatments at the $p < 0.05$ level according to Duncan’s multiple range test.

3.4. Analysis of Related Gene Expression Patterns during the Regeneration of *F. mandshurica*

The expression levels of the *FmYUC1*, *FmWOX4*, *FmWOX5*, *FmWUS*, *FmARR2*, *FmREV*, *FmPLT2* and *FmPLT3* genes increased significantly in the early stage of adventitious bud induction culture (Figure 8). The expression of *FmYUC1*, *FmPLT2* and *FmPLT3* cultivated for 1 d was 48.16, 5.93 and 1.57 times the 0 d cultivated level. At the stage of callus formation at 1–15 d, the *ARR2* gene had a higher expression level than the control. The expression levels of *FmYUC1*, *FmWOX4*, *FmWOX5*, *FmWUS*, *FmARR2*, *FmREV*, *FmZPR3*, *FmPLT2* and *FmPLT3* were markedly increased at 1 d. The expression of *FmWUS*, *FmWOX4*, *FmWOX5* and *FmARR2* cultivated for 1 d was 1.67, 5.69, 2.28 and 6.56 times the 0 d cultivated level. The expression of *FmWIND1*, *FmYUC1*, *FmWOX5*, *FmREV* and *FmPLT3* was greater than the control group after callus formation, and the adventitious buds were formed at 30 d. The expression of the HD-Zip III transcription factor *FmREV* cultivated for 1 d was 10 times the 0 d cultivated level. The expression of *FmWIND1*, *FmYUC1*, *FmWOX5*, *FmREV* and *FmPLT3* was higher than that of the control, indicating that these genes are associated with the process of adventitious buds.




Figure 8. Gene expression of adventitious buds in the stem segments of *F. mandshurica*. (a–i) 0–15 d was used to induce callus in the dark culture stage. After 15–30 d, it was transferred to the light culture stage to induce callus regeneration, and the shoot regeneration time was 30 d. Different letters above bars within statistically significant differences between different times of the treatments at the $p < 0.05$ level according to Duncan's multiple range test.

3.5. Establishment of a Callus Regeneration System from Stem Segments of *F. mandshurica*

Stem tissue of *F. mandshurica* was used as the explant (Figure 9a). Stem tissue was cultured for 20 d until it formed a callus (Figure 9b). The induction rate of the callus was 99.15%. Stem tissue

was cultured for 30 d until it formed adventitious buds (Figure 9c), and the induction rate of the adventitious bud was 51.11%. After transferring to the secondary proliferation medium (Figure 9d–f), the shoots were transferred to the rooting medium for three weeks. The rooting rate was 80.54% (Figure 9g). After alternating light and dark cultivation, the survival rate reached 90% (Figure 9h). Finally, the viable *F. mandshurica* was transferred into the pot to continue cultivation (Figure 9i), and the in vitro regeneration technique of *F. mandshurica* was obtained.




Figure 9. Preliminary establishment of the callus regeneration system of stems of *F. mandshurica*. (a) Inoculation of stem explants (scale 5 mm); (b) Callus formed at 20 d; (c) at 30 d, young adventitious buds form; (d) at 40 d, clustered buds formed; (e) Adventitious bud elongation; (f) Sub-generation expansion; (g) Rooting (scale 5 mm); (h) Planting (scale 1 cm); (i) Sterilization of sterile *F. mandshurica* seedlings to survive (scale 3 cm).

4. Discussion

4.1. Effect of the Exogenous Hormones TDZ and 6-BA on Induction of Adventitious Buds of the Callus

The experimental results showed that, at the same TDZ concentration, when 0–8 mg/L 6-BA was added, the higher the concentration of 6-BA, the looser the callus. However, when the adventitious buds induced by cotyledons were used as explants, a higher concentration of 6-BA did not contribute to the formation of a cotyledon callus [23]. The texture of the callus formed by the low-concentration 6-BA-induced material was dense and opaque green, and the texture surface of the callus formed by the high-concentration 6-BA-induced material was relatively loose and emerald green (Figure 1a,b).

The callus formed by the hypocotyl was rigid (Figure 1d). At the same 6-BA concentration, with an increase in TDZ concentration, the callus induction rate increased, indicating that TDZ plays a key role in hypocotyl callus induction. This is consistent with the results of Zeng et al., who found that the frequency of bud organ development of explants cultured with TDZ (71.16%) is far higher than that of explants cultured without TDZ (30.42%) [24]. Internal callus tissue began to age and darken after 30 days of culture, and did not contribute to the formation of adventitious buds. In addition, callus induced by using roots as explants was milky white in the lower culture and yellow-green in the upper culture. After subculture, the internal callus was obviously browned, while the external callus could not produce a primordium structure and adventitious buds (Figure 1c). The callus formed by cotyledons has the characteristics of harder and slower growth, and the callus formed by roots has the characteristics of loose and fast growth. This may be why the cotyledons and roots failed to induce adventitious buds.

4.2. Wound Treatment Is Beneficial for Callus Re-Differentiation

An important aspect of organ regeneration is the restoration of cell proliferation and regenerative capacity at the wound site, which results in the formation of callus and subsequent establishment of shoot or root apical meristems [25]. In this study, the research indicates that the callus formed at the wound site of the stem segment differed from the callus formed by the epidermis of the stem segment and was more conducive to the formation and differentiation of adventitious buds (Figure 2). Studies have found that Arabidopsis leaves may produce at least two wound signals after trauma, triggering short- and long-term wound signals, respectively. Short-term wound signals are mainly involved in controlling the fate transition of auxin and regenerative cells, while long-term wound signals primarily regulate the cellular environment at the wound site and maintain auxin levels to form regenerative cells [15]. The expression levels of *FmYUC1*, *FmWOX4*, *FmWOX5*, *FmWUS*, *FmARR2*, *FmREV*, *FmPLT2* and *FmPLT3* increased to different degrees after initial access to callus induction medium, the auxin synthesis gene *YUC1* and the auxin pathway genes. *FmPLT2* and *FmPLT3* had an expression level of 48.16, 5.93, and 1.57 times the control. The expression levels of the somatic embryogenesis and regeneration genes *FmWUS*, *FmWOX4* and *FmWOX5* were 1.67, 5.69 and 2.28 times the control. The expression level of the cytokinin pathway gene b-type ARRs transcription factor *FmARR2* in Arabidopsis shoot apical meristems was 6.56 times that of the control [26]. The expression level of the HD-ZIP III transcription factor *FmREV* was 0.4 and 10.84 times that of the control. At 1–15 d, at the callus formation stage, the expression level of the *FmARR2* gene at 15 d was higher than that at 0 d, indicating that the gene is involved in the regulation of callus formation, and the gene responds to cytokinin signals [27]. The results showed that the formation of callus may be regulated by cytokinin signal transduction. The expression of *FmWIND1*, *FmYUC1*, *FmWOX5*, *FmREV* and *FmPLT3* was up-regulated compared with the control during the formation of adventitious buds from 15–30 days. These genes may be involved in the formation of buds. *FmREV*, *FmPLT3* and *FmYUC1* belong to the process of bud regeneration. In the regulation of auxin synthesis genes [28], *FmWIND1* is involved in the formation of cytokinin pathways and cells in the programming process [29,30], and *FmWOX5* can maintain callus embryogenicity [31]. The occurrence of adventitious buds may be affected by the *FmWIND1* gene induced by the wound signal, which changes the fate of the cells [32–34].

4.3. 5-aza Promotes Callus Re-Differentiation

In the plant genome, DNA methylation is mostly 5-mC and 7-mG. During plant tissue culture, different species, periods, genotypes, hormones, and DNA methylation inhibitors can make DNA methylation levels different [35–38]. Low DNA methylation levels favor Arabidopsis root regeneration [20]. During the process of callus formation and regeneration of sprouts, DNA methylation presents a variety of trends [39–41]. Dedifferentiated maize embryos have higher levels of methylation than normal embryos, and DNA methylation levels increase during callus induction [42]. The molecular mechanism involved in DNA methylation and regeneration is primarily focused on the study of

FmWUS, but unfortunately, in this study, there was no change in the expression of the *FmWUS* gene, which may be that the cell range of *FmWUS* expression is small, and the expression time is short. The number of adventitious buds of *F. mandshurica* is small, the expression level of multi-cells may be averaged, and the expression of *FmWUS* is masked. However, the experimental results showed that the expression levels of *FmREV*, *FmYUC1* and *FmWOX4* increased after 5-aza was added, which indicated that 5-aza could promote the expression of these three genes and lead to an increased rate of adventitious bud differentiation.

5. Conclusions

The purpose of this study was to establish and optimize efficient and effective callus induction and an adventitious bud regeneration system. By comparing the differences between callus induction and adventitious bud differentiation of four kinds of explants and combining epigenetic regulation analysis, the regeneration technology of explant callus from stem segments of *F. mandshurica* were established and optimized. The changes in 5-aza and HDAC enzyme activity in callus formation and adventitious bud differentiation were explored. The level of DNA methylation in the callus formation stage of *F. mandshurica* decreased, and the adventitious bud formation stage increased again. HDAC enzyme activity was significantly reduced during the induction of adventitious buds, and TSA inhibited the activity of HDAC enzymes to promote adventitious bud regeneration. Analysis of the endogenous hormone content and related gene expression during adventitious bud differentiation showed that exogenous 6-BA and TDZ promote the differentiation of adventitious buds in governing the levels of endogenous IAA and ZR hormones. *FmWIND1*, *FmWOX5* and *FmPLT3* advance the formation of adventitious buds.

Author Contributions: L.Y. and Y.Z. designed the experiment and performed technical guidance. L.Y. wrote the manuscript and provided assistance in the data analysis. All authors participated in the experimental process. All authors have read and agreed to the published version of the manuscript.

Funding: This study was funded supported by the Applied Technology Research and Developmental Plan Program of Heilongjiang Province (GA19B201).

Conflicts of Interest: The authors declare no conflict of interest.

References

- Gu, J.C.; Wei, X.; Wang, J.; Dong, X.Y.; Wang, Z.Q. Marked differences in standing biomass, length density, anatomy and physiological activity between white and brown roots in *Fraxinus mandshurica* rupr. plantation. *Plant Soil* **2015**, *392*, 267–277. [\[CrossRef\]](#)
- Wang, J.S.; Wu, L.; Zhao, X.H.; Fan, J.; Zhang, C.Y.; von Gadow, K. Influence of ground flora on *Fraxinus mandshurica* seedling growth on abandoned land and beneath forest canopy. *Eur. J. For. Res.* **2013**, *132*, 313–324. [\[CrossRef\]](#)
- Chen, X.H.; Wang, H.T.; Li, X.Y.; Ma, K.; Zhan, Y.G.; Zeng, F.S. Molecular cloning and functional analysis of 4-coumarate:Coa ligase 4(4cl-like 1) from *Fraxinus mandshurica* and its role in abiotic stress tolerance and cell wall synthesis. *BMC Plant Biol.* **2019**, *19*, 201. [\[CrossRef\]](#) [\[PubMed\]](#)
- Liu, W.G.; Fan, X.H.; Wang, J.S.; Zhang, C.Y.; Lu, W.M.; Gadow, K.V. Spectral reflectance response of *Fraxinus mandshurica* leaves to above- and belowground competition. *Int. J. Remote Sens.* **2012**, *33*, 5072–5086. [\[CrossRef\]](#)
- Sun, X.D.; Wang, W.N.; Razaq, M.; Sun, H.L. Effects of stand density on soil organic carbon storage in the top and deep soil layers of *Fraxinus mandshurica* plantations. *Austrian J. For. Sci.* **2019**, *136*, 27–44.
- He, L.M.; Zeng, F.S.; He, Z.L.; Li, Y.H.; Li, S.J.; Zhan, Y.G. Interspecific hybridizations of *Fraxinus* l. (*F. mandshurica* × *f. Americana* and *F. mandshurica* × *f. Velutina*) and heterosis analysis and selection of *f1* progenies. *Can. J. For. Res.* **2019**, *49*, 1265–1276. [\[CrossRef\]](#)
- Kong, D.M.; Preece, J.E.; Shen, H.L. Somatic embryogenesis in immature cotyledons of manchurian ash (*Fraxinus mandshurica* rupr.). *Plant Cell Tissue Organ. Cult.* **2012**, *108*, 485–492. [\[CrossRef\]](#)

8. Stener, L.G. Genetic evaluation of damage caused by ash dieback with emphasis on selection stability over time. *For. Ecol. Manag.* **2018**, *409*, 584–592. [[CrossRef](#)]
9. Inoue, T.; Okane, I.; Ishiga, Y.; Degawa, Y.; Hosoya, T.; Yamaoka, Y. The life cycle of *hymenoscyphus fraxineus* on manchurian ash, *fraxinus mandshurica*, in Japan. *Mycoscience* **2019**, *60*, 89–94. [[CrossRef](#)]
10. Du, X.Z.; Wei, X. Definition of fine roots on the basis of the root anatomy, diameter, and branch orders of one-year old *fraxinus mandshurica* seedlings. *J. For. Res.* **2018**, *29*, 1321–1327. [[CrossRef](#)]
11. Lee, J.H.; Pijut, P.M. Adventitious shoot regeneration from in vitro leaf explants of *fraxinus nigra*. *Plant Cell Tissue Organ. Cult.* **2017**, *130*, 335–343. [[CrossRef](#)]
12. Aida, M.; Beis, D.; Heidstra, R.; Willemsen, V.; Blilou, I.; Galinha, C.; Nussaume, L.; Noh, Y.S.; Amasino, R.; Scheres, B. The *plethora* genes mediate patterning of the arabidopsis root stem cell niche. *Cell* **2004**, *119*, 109–120. [[CrossRef](#)] [[PubMed](#)]
13. Iwase, A.; Mitsuda, N.; Koyama, T.; Hiratsu, K.; Kojima, M.; Arai, T.; Inoue, Y.; Seki, M.; Sakakibara, H.; Sugimoto, K.; et al. The *ap2/erf* transcription factor *wind1* controls cell dedifferentiation in arabidopsis. *Curr. Biol.* **2011**, *21*, 508–514. [[CrossRef](#)] [[PubMed](#)]
14. Zhou, C.; Guo, J.S.; Feng, Z.H.; Cui, X.H.; Zhu, J. Molecular characterization of a novel *ap2* transcription factor *thwind1-1* from *thellungiella halophila*. *Plant Cell Tissue Organ. Cult.* **2012**, *110*, 423–433. [[CrossRef](#)]
15. Chen, L.Q.; Sun, B.B.; Xu, L.; Liu, W. Wound signaling: The missing link in plant regeneration. *Plant Signal. Behav.* **2016**, *11*, 4273–4284. [[CrossRef](#)] [[PubMed](#)]
16. Xin, W.; Wang, Z.C.; Liang, Y.; Wang, Y.H.; Hu, Y.X. Dynamic expression reveals a two-step patterning of *wus* and *clv3* during axillary shoot meristem formation in arabidopsis. *J. Plant Physiol.* **2017**, *214*, 1–6. [[CrossRef](#)]
17. Kareem, A.; Durgaprasad, K.; Sugimoto, K.; Du, Y.J.; Pulianmackal, A.J.; Trivedi, Z.B.; Abhayadev, P.V.; Pinon, V.; Meyerowitz, E.M.; Scheres, B.; et al. *Plethora* genes control regeneration by a two-step mechanism. *Curr. Biol.* **2015**, *25*, 1017–1030. [[CrossRef](#)]
18. Wojcikowska, B.; Botor, M.; Moronczyk, J.; Wojcik, A.M.; Nodzynski, T.; Karcz, J.; Gaj, M.D. Trichostatin a triggers an embryogenic transition in arabidopsis explants via an auxin-related pathway. *Front. Plant Sci.* **2018**, *9*, 1353. [[CrossRef](#)]
19. Zhao, Q.; Du, Y.L.; Wang, H.T.; Rogers, H.J.; Yu, C.M.; Liu, W.; Zhao, M.Z.; Xie, F.T. 5-azacytidine promotes shoot regeneration during agrobacterium-mediated soybean transformation. *Plant Physiol. Biochem.* **2019**, *141*, 40–50. [[CrossRef](#)]
20. Shemer, O.; Landau, U.; Candela, H.; Zemach, A.; Williams, L.E. Competency for shoot regeneration from arabidopsis root explants is regulated by DNA methylation. *Plant Sci.* **2015**, *238*, 251–261. [[CrossRef](#)]
21. Liang, N.S.; Zhan, Y.G.; Yu, L.; Wang, Z.Q.; Zeng, F.S. Characteristics and expression analysis of *fmtcp15* under abiotic stresses and hormones and interact with *della* protein in *fraxinus mandshurica* rupr. *Forests* **2019**, *10*, 343. [[CrossRef](#)]
22. Liang, N.S.; Yu, L.; Liu, C.H.; Wang, Z.Q.; Zhao, X.T.; Zhan, Y.G. Molecular cloning and expression under abiotic stresses and hormones of the ethylene response factor *vii* gene *fmrp2.12* from *fraxinus mandshurica*. *J. For. Res.* **2019**, *30*, 1289–1300. [[CrossRef](#)]
23. Yang, X.; Yang, X.; Guo, T.; Gao, K.; Zhao, T.; Chen, Z.; An, X. High-efficiency somatic embryogenesis from seedlings of *koelreuteria paniculata* laxm. *Forests* **2018**, *9*, 769. [[CrossRef](#)]
24. Zeng, Q.Q.; Han, Z.Q.; Kang, X.Y. Adventitious shoot regeneration from leaf, petiole and root explants in triploid (*populus alba*x*p.glandulosa*)x*p.tomentosa*. *Plant Cell Tissue Organ. Cult.* **2019**, *138*, 121–130. [[CrossRef](#)]
25. Ikeuchi, M.; Shibata, M.; Rymen, B.; Iwase, A.; Bagman, A.M.; Watt, L.; Coleman, D.; Favero, D.S.; Takahashi, T.; Ahnert, S.E.; et al. A gene regulatory network for cellular reprogramming in plant regeneration. *Plant Cell Physiol.* **2018**, *59*, 770–782. [[CrossRef](#)]
26. Kim, Y.S.; Kim, S.G.; Lee, M.; Lee, I.; Park, H.Y.; Seo, P.J.; Jung, J.H.; Kwon, E.J.; Suh, S.W.; Paek, K.H.; et al. *Hd-zip iii* activity is modulated by competitive inhibitors via a feedback loop in arabidopsis shoot apical meristem development. *Plant Cell* **2008**, *20*, 920–933. [[CrossRef](#)]
27. Shull, T.E.; Kurepa, J.; Smalle, J.A. Cytokinin signaling promotes differential stability of type-b *arrs*. *Plant Signal. Behav.* **2016**, *11*, e1169354. [[CrossRef](#)]
28. Liu, Y.; Zhang, H.L.; Guo, H.R.; Xie, L.; Zeng, R.Z.; Zhang, X.Q.; Zhang, Z.S. Transcriptomic and hormonal analyses reveal that *yuc*-mediated auxin biogenesis is involved in shoot regeneration from rhizome in *cymbidium*. *Front. Plant Sci.* **2017**, *8*. [[CrossRef](#)]

29. Ikeda, M.; Ohme-Takagi, M. Tcps, wuss, and winds: Families of transcription factors that regulate shoot meristem formation, stem cell maintenance, and somatic cell differentiation. *Front. Plant Sci.* **2014**, *5*. [[CrossRef](#)]
30. Iwase, A.; Mita, K.; Favero, D.S.; Mitsuda, N.; Sasaki, R.; Kobayshi, M.; Takebayashi, Y.; Kojima, M.; Kusano, M.; Oikawa, A.; et al. Wind1 induces dynamic metabolomic reprogramming during regeneration in brassica napus. *Dev. Biol.* **2018**, *442*, 40–52. [[CrossRef](#)]
31. Orłowska, A.; Kepczynska, E. Identification of polycomb repressive complex1, trithorax group genes and their simultaneous expression with wuschel, wuschel-related homeobox5 and shoot meristemless during the induction phase of somatic embryogenesis in medicago truncatula gaertn. *Plant Cell Tissue Organ. Cult.* **2018**, *134*, 345–356. [[CrossRef](#)]
32. Iwase, A.; Mitsuda, N.; Ikeuchi, M.; Ohnuma, M.; Koizuka, C.; Kawamoto, K.; Imamura, J.; Ezura, H.; Sugimoto, K. Arabidopsis wind1 induces callus formation in rapeseed, tomato, and tobacco. *Plant Signal. Behav.* **2013**, *8*, e27432. [[CrossRef](#)] [[PubMed](#)]
33. Iwase, A.; Mita, K.; Nonaka, S.; Ikeuchi, M.; Koizuka, C.; Ohnuma, M.; Ezura, H.; Imamura, J.; Sugimoto, K. Wind1-based acquisition of regeneration competency in arabidopsis and rapeseed. *J. Plant Res.* **2015**, *128*, 389–397. [[CrossRef](#)] [[PubMed](#)]
34. Iwase, A.; Harashima, H.; Ikeuchi, M.; Rymen, B.; Ohnuma, M.; Komaki, S.; Morohashi, K.; Kurata, T.; Nakata, M.; Ohme-Takagi, M.; et al. Wind1 promotes shoot regeneration through transcriptional activation of enhancer of shoot regeneration1 in arabidopsis. *Plant Cell* **2017**, *29*, 54–69. [[CrossRef](#)]
35. Salaj, T.; Klubicova, K.; Matusova, R.; Salaj, J. Somatic embryogenesis in selected conifer trees pinus nigra arn. And abies hybrids. *Front. Plant Sci.* **2019**, *10*. [[CrossRef](#)]
36. Ree, J.F.; Polesi, L.G.; Back, F.; Bertolazi, A.A.; Silveira, V.; Guerra, M.P. Aging peach palm (bactris gasipaes kunth) cultures lose embryogenic potential and metabolic cellular function due to continuous culture in hypoxic environments. *Plant Cell Tissue Organ. Cult.* **2019**, *140*, 49–67. [[CrossRef](#)]
37. Muniandy, K.; Tan, M.H.; Song, B.K.; Ayub, Q.; Rahman, S. Comparative sequence and methylation analysis of chloroplast and amyloplast genomes from rice. *Plant Mol. Biol.* **2019**, *100*, 33–46. [[CrossRef](#)]
38. Li, J.Y.; Wang, M.J.; Li, Y.J.; Zhang, Q.H.; Lindsey, K.; Daniell, H.; Jin, S.X.; Zhang, X.L. Multi-omics analyses reveal epigenomics basis for cotton somatic embryogenesis through successive regeneration acclimation process. *Plant Biotechnol. J.* **2019**, *17*, 435–450. [[CrossRef](#)]
39. Hu, L.J.; Li, N.; Zhang, Z.B.; Meng, X.C.; Dong, Q.L.; Xu, C.M.; Gong, L.; Liu, B. Cg hypomethylation leads to complex changes in DNA methylation and transpositional burst of diverse transposable elements in callus cultures of rice. *Plant J.* **2019**, *101*, 188–203. [[CrossRef](#)]
40. Gulyas, A.; Dobranszki, J.; Kiss, E.; da Silva, J.A.T.; Posta, K.; Hidvegi, N. Changes in DNA methylation pattern of apple long-term in vitro shoot culture and acclimatized plants. *J. Plant Physiol.* **2019**, *239*, 18–27. [[CrossRef](#)]
41. Osorio-Montalvo, P.; Saenz-Carbonell, L.; De-la-Pena, C. 5-azacytidine: A promoter of epigenetic changes in the quest to improve plant somatic embryogenesis. *Int. J. Mol. Sci.* **2018**, *19*, 3182. [[CrossRef](#)] [[PubMed](#)]
42. Liu, H.J.; Ma, L.L.; Yang, X.R.; Zhang, L.; Zeng, X.; Xie, S.P.; Peng, H.W.; Gao, S.B.; Lin, H.J.; Pan, G.T.; et al. Integrative analysis of DNA methylation, mrnas, and small rnas during maize embryo dedifferentiation. *BMC Plant Biol.* **2017**, *17*. [[CrossRef](#)] [[PubMed](#)]

