

## Anti-ADAMTS13 antibody

Cat. No.	ml161236
Package	25 μl/100 μl/200 μl
Storage	-20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-ADAMTS13 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human ADAMTS13
Reactivity	Human
Content	0.6 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	ADAMTS13

rarget information	
Symbol	ADAMTS13
Full name	ADAM metallopeptidase with thrombospondin type 1 motif, 13
Synonyms	VWFCP; C9orf8; vWF-CP; ADAM-TS13; ADAMTS-13
Swissprot	Q76LX8
Target Background	

,00<sup>c</sup>


## **Target Background**

This gene encodes a member of a family of proteins containing several distinct regions, including a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. The enzyme encoded by this gene specifically cleaves von Willebrand Factor (vWF) Defects in this gene are associated with thrombotic thrombocytopenic purpura. Alternative splicing results in multiple transcript variants.



## Applications Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human liver cancer Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human liver cancer tissue using ml161236(ADAMTS13 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 2000-5000

Predicted cell location: Cytoplasm Positive control: Human colorectal cancer Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml161236(ADAMTS13 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696
- 邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn