

订购热线: 4008-898-798

# Anti-AKT2 antibody

**Cat. No.** ml161251

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-AKT2 rabbit polyclonal antibody

**Applications** ELISA, IHC

**Immunogen** Synthetic peptide of human AKT2

Reactivity Human, Mouse, Rat

Content 1 mg/ml Host species Rabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol AKT2

Full name v-akt murine thymoma viral oncogene homolog 2
Synonyms PKBB; PRKBB; HIHGHH; PKBBETA; RAC-BETA

Swissprot P31751

### **Target Background**

This gene is a putative oncogene encoding a protein belonging to a subfamily of serine/threonine kinases containing SH2-like (Src homology 2-like) domains. The gene was shown to be amplified and overexpressed in 2 of 8 ovarian carcinoma cell lines and 2 of 15 primary ovarian tumors. Overexpression contributes to the malignant phenotype of a subset of human ductal pancreatic cancers. The encoded protein is a general protein kinase capable of phophorylating several known proteins.



订购热线: 4008-898-798


# **Applications**


## **Immunohistochemistry**

Predicted cell location: Cytoplasm

Positive control: Human esophagus cancer

Recommended dilution: 30-150





The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml161251(AKT2 Antibody) at dilution 1/40, on the right is treated with synthetic peptide. (Original magnification: ×200)

#### **ELISA**

Recommended dilution: 2000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn