

订购热线: 4008-898-798

# Anti-VGLL3 antibody

**Cat. No.** ml163165

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-VGLL3 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human VGLL3

ReactivityHuman, MouseContent0.6 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol VGLL3

Full name vestigial like family member 3

Synonyms VGL3; VGL-3 Swissprot A8MV65

## **Target Background**

May act as a specific coactivator for the mammalian TEFs.



订购热线: 4008-898-798


### **Applications**


#### **Immunohistochemistry**

Predicted cell location: Nucleus

Positive control: Human colorectal cancer

Recommended dilution: 20-100






The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using ml163165(VGLL3 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

Predicted cell location: Nucleus

Positive control: Human gastric cancer

Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using ml163165(VGLL3 Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

### **ELISA**

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn