

Anti-UBR4 antibody

Cat. No.	ml163223
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-UBR4 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Synthetic peptide of human UBR4
Reactivity	Human, Mouse, Rat
Content	1.1 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	UBR4
Full name	ubiquitin protein ligase E3 component n-recognin 4
Synonyms	p600; ZUBR1; RBAF600
Swissprot	Q5T4S7

Target Background


The protein encoded by this gene is an E3 ubiquitin-protein ligase that interacts with the retinoblastoma-associated protein in the nucleus and with calcium-bound calmodulin in the cytoplasm. The encoded protein appears to be a cytoskeletal component in the cytoplasm and part of the chromatin scaffold in the nucleus. In addition, this protein is a target of the human papillomavirus type 16 E7 oncoprotein.

订购热线: 4008-898-798


Applications

Immunohistochemistry


Predicted cell location: Cytoplasm
Positive control: Human brain
Recommended dilution: 40-200




The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml163223(UBR4 Antibody) at dilution 1/65, on the right is treated with synthetic peptide. (Original magnification: ×200)



Predicted cell location: Cytoplasm
Positive control: Human lung cancer
Recommended dilution: 40-200



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml163223(UBR4 Antibody) at dilution 1/65, on the right is treated with synthetic peptide. (Original magnification: ×200)



ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn