

订购热线: 4008-898-798

Anti-DHRS4 antibody

Cat. No. ml163712

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-DHRS4 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human DHRS4

ReactivityHumanContent0.48 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol DHRS4

Full name dehydrogenase/reductase 4

Synonyms CR; NRDR; PHCR; PSCD; SDR-SRL; SDR25C1; SDR25C2; SCAD-SRL

Swissprot Q9BTZ2

Target Background

Reduces all-trans-retinal and 9-cis retinal. Can also catalyze the oxidation of all-trans-retinol with NADP as co-factor, but with much lower efficiency. Reduces alkyl phenyl ketones and alpha-dicarbonyl compounds with aromatic rings, such as pyrimidine-4-aldehyde, 3-benzoylpyridine, 4-benzoylpyridine, menadione and 4-hexanoylpyridine. Has no activity towards aliphatic aldehydes and ketones (By similarity).




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm Positive control: Human ovarian cancer

Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml163712(DHRS4 Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn