

订购热线: 4008-898-798

Anti-POLR2L antibody

Cat. No. ml164037

Package 25 μ l/100 μ l/200 μ l

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-POLR2L rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human POLR2L

ReactivityHuman, MouseContent1.32 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol POLR2L

Full name RNA polymerase II subunit L

Synonyms RBP10; RPB10; RPABC5; RPB7.6; hRPB7.6; RPB10beta

Swissprot P62875

Target Background

This gene encodes a subunit of RNA polymerase II, the polymerase responsible for synthesizing messenger RNA in eukaryotes. The product of this gene contains four conserved cysteines characteristic of an atypical zinc-binding domain. Like its counterpart in yeast, this subunit may be shared by the other two DNA-directed RNA polymerases.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm or Nucleus Positive control: Human cervical cancer

Recommended dilution: 40-200





The image on the left is immunohistochemistry of paraffin-embedded Human cervical cancer tissue using ml164037(POLR2L Antibody) at dilution 1/45, on the right is treated with synthetic peptide. (Original

magnification: ×200)

ELISA


Recommended dilution: 5000-10000


联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn

Predicted cell location: Cytoplasm or Nucleus Positive control: Human thyroid cancer Recommended dilution: 40-200





The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml164037(POLR2L Antibody) at dilution 1/45, on the right is treated with synthetic peptide. (Original magnification: ×200)