

订购热线: 4008-898-798

Anti-GPRC5C antibody

Cat. No. ml164117

Package 25 μl/100 μl/200 μl

Storage -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview

Description Anti-GPRC5C rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Synthetic peptide of human GPRC5C

ReactivityHuman, MouseContent0.78 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

Target information

Symbol GPRC5C

Full name G protein-coupled receptor class C group 5 member C

Synonyms RAIG3; RAIG-3

Swissprot Q9NQ84

Target Background

The protein encoded by this gene is a member of the type 3 G protein-coupled receptor family. Members of this superfamily are characterized by a signature 7-transmembrane domain motif. The specific function of this protein is unknown; however, this protein may mediate the cellular effects of retinoic acid on the G protein signal transduction cascade. Two transcript variants encoding different isoforms have been found for this gene.




订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm and Cell membrane

Positive control: Human tonsil Recommended dilution: 20-100





The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml164117(GPRC5C Antibody) at dilution 1/30, on the right is treated with synthetic peptide. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net 网址: www.mlbio.cn