

## Anti-LGALSL antibody

<b>Cat. No.</b>	ml164131
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-LGALSL rabbit polyclonal antibody
<b>Applications</b>	ELISA, IHC
<b>Immunogen</b>	Synthetic peptide of human LGALSL
<b>Reactivity</b>	Human, Mouse
<b>Content</b>	1.68 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	LGALSL
<b>Full name</b>	galectin like
<b>Synonyms</b>	GRP; HSPC159
<b>Swissprot</b>	Q3ZCW2

### Target Background



Does not bind lactose, and may not bind carbohydrates.

订购热线: 4008-898-798

### Applications



#### Immunohistochemistry

Predicted cell location: Nucleus and Cytoplasm  
Positive control: Human ovarian cancer  
Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml164131(LGALS1 Antibody) at dilution 1/70, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )

Predicted cell location: Nucleus and Cytoplasm  
Positive control: Human thyroid cancer  
Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml164131(LGALS1 Antibody) at dilution 1/70, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )

#### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)