

## Anti-HPDL antibody

<b>Cat. No.</b>	ml164214
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-HPDL rabbit polyclonal antibody
<b>Applications</b>	ELISA, IHC
<b>Immunogen</b>	Synthetic peptide of human HPDL
<b>Reactivity</b>	Human
<b>Content</b>	0.6 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	HPDL
<b>Full name</b>	4-hydroxyphenylpyruvate dioxygenase like
<b>Synonyms</b>	GLOXD1; 4-HPPD-L
<b>Swissprot</b>	Q96IR7

### Target Background


May have dioxygenase activity.

订购热线: 4008-898-798


### Applications

#### Immunohistochemistry


Predicted cell location: Cytoplasm  
Positive control: Human tonsil  
Recommended dilution: 25-100




The image on the left is immunohistochemistry of paraffin-embedded Human tonsil tissue using ml164214(HPDL Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )



Predicted cell location: Cytoplasm  
Positive control: Human brain  
Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml164214(HPDL Antibody) at dilution 1/20, on the right is treated with synthetic peptide. (Original magnification:  $\times 200$ )



#### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)