

Anti-OASL antibody

Cat. No.	ml221065
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-OASL rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human OASL
Reactivity	Human
Content	0.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	OASL
Full name	2'-5'-oligoadenylate synthetase-like
Synonyms	OASLd, TRIP14, TRIP-14, p59OASL, p59 OASL, p59-OASL
Swissprot	Q15646

Target Background

OASL encodes an interferon-induced protein that is a member of the 2-5A synthetase family. 2-prime,5-prime oligoadenylates (2-5As) bind to and activate RNase L, resulting in general RNA degradation and consequent inhibition of protein synthesis. 2-5As are produced by a well-conserved family of interferon-induced enzymes, the 2-5A synthetases or OASs

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human colon cancer

Recommended dilution: 25-100



The image on the left is immunohistochemistry of paraffin-embedded Human colon cancer tissue using ml221065(OASL Antibody) at dilution 1/35, on the right is treated with fusion protein. (Original magnification: $\times 200$)

ELISA

Recommended dilution: 2000-5000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn