

订购热线: 4008-898-798

# Anti-RCHY1 antibody

**Cat. No.** ml221133

**Package** 25 μl/100 μl/200 μl

**Storage** -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

**Product overview** 

**Description** Anti-RCHY1 rabbit polyclonal antibody

Applications ELISA, IHC

Immunogen Fusion protein of human RCHY1

ReactivityHuman, MouseContent0.2 mg/mlHost speciesRabbit

Ig classImmunogen-specific rabbit IgGPurificationAntigen affinity purification

**Target information** 

Symbol RCHY1

Full name ring finger and CHY zinc finger domain containing 1, E3 ubiquitin protein ligase

Synonyms ARNIP, CHIMP, PIRH2, PIRH2E, PIRH2F, RNF199, ZNF363, hARNIP, PRO1996

Swissprot Q96PM5

### **Target Background**

The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including TP53, HDAC1 and CDKN1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.




订购热线: 4008-898-798


### **Applications**

## **Immunohistochemistry**

Predicted cell location: Cytoplasm Positive control: Human ovarian cancer

Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human ovarian cancer tissue using ml221133(RCHY1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original

magnification: ×200)

#### **ELISA**


Recommended dilution: 1000-5000


联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn

Predicted cell location: Cytoplasm Positive control: Human lung cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml221133(RCHY1 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)