

Anti-ACOT7 antibody

Cat. No.	ml221469
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-ACOT7 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human ACOT7
Reactivity	Human
Content	0.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	ACOT7
Full name	acyl-CoA thioesterase 7
Synonyms	ACT; ACH1; BACH; LACH; LACH1; hBACH; CTE-II
Swissprot	O00154

Target Background

This gene encodes a member of the acyl coenzyme family. The encoded protein hydrolyzes the CoA thioester of palmitoyl-CoA and other long-chain fatty acids. Decreased expression of this gene may be associated with mesial temporal lobe epilepsy. Alternatively spliced transcript variants encoding distinct isoforms with different subcellular locations have been characterized.

订购热线: 4008-898-798


Applications

Immunohistochemistry

Predicted cell location: Cytoplasm

Positive control: Human brain

Recommended dilution: 20-100



The image is immunohistochemistry of paraffin-embedded Human brain tissue using ml221469(ACOT7 Antibody) at dilution 1/25.
(Original magnification: ×200)

ELISA

Recommended dilution: 1000-2000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn