

## Anti-COMMD6 antibody

 Cat. No.
 ml123771

 Package
 25 μl/100 μl/200 μl

 Storage
 -20°C, pH7.4 PBS, 0.05% NaN3, 40% Glycerol

Product overview	
Description	Anti-COMMD6 rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Full length fusion protein
Reactivity	Human, Mouse
Content	0.9 mg/ml
Host species	Rabbit
lg class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification
Target information	
Symbol	COMMD6

COMMD6 COMM domain containing 6 Acrg Q7Z4G1

## **Target Background**

Full name

Synonyms

Swissprot


COMMD6 belongs to a family of NF-kappa-B (see RELA; MIM 164014)-inhibiting proteins characterized by the presence of a COMM domain (see COMMD1; MIM 607238) (de Bie et al., 2006 [PubMed 16573520])




订购热线: 4008-898-798

## Applications Immunohistochemistry

Predicted cell location: Nucleus Positive control: Human prostate cancer Recommended dilution: 25-100







The image on the left is immunohistochemistry of paraffin-embedded Human prostate cancer tissue using ml123771(COMMD6 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

## ELISA

Recommended dilution: 5000-10000

Predicted cell location: Nucleus Positive control: Human esophagus cancer Recommended dilution: 25-100





The image on the left is immunohistochemistry of paraffin-embedded Human esophagus cancer tissue using ml123771(COMMD6 Antibody) at dilution 1/30, on the right is treated with fusion protein. (Original magnification: ×200)

- 联系电话: 4008-898-798, 021-61725725
- 联系QQ: 2881505695,2881505696
- 邮箱: mlbio\_cn@yeah.net 网址: www.mlbio.cn