

Anti-VILL antibody

Cat. No.	ml124109
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-VILL rabbit polyclonal antibody
Applications	ELISA, IHC
Immunogen	Fusion protein of human VILL
Reactivity	Human
Content	1.2 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	VILL
Full name	villin-like
Synonyms	
Swissprot	O15195

Target Background



The protein encoded by this gene belongs to the villin/gelsolin family. It contains 6 gelsolin-like repeats and a headpiece domain. It may play a role in actin-bundling.

订购热线: 4008-898-798

Applications



Immunohistochemistry

Predicted cell location: Cytoplasm
Positive control: Human brain
Recommended dilution: 50-350



The image on the left is immunohistochemistry of paraffin-embedded Human brain tissue using ml124109(VILL Antibody) at dilution 1/65, on the right is treated with fusion protein. (Original magnification: ×200)

Predicted cell location: Cytoplasm
Positive control: Human lung cancer
Recommended dilution: 50-350



The image on the left is immunohistochemistry of paraffin-embedded Human lung cancer tissue using ml124109(VILL Antibody) at dilution 1/65, on the right is treated with fusion protein. (Original magnification: ×200)

ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn