

## Anti-BRD9 antibody

**Cat. No.** ml124272  
**Package** 25 µl/100 µl/200 µl  
**Storage** -20°C, pH7.4 PBS, 0.05% NaN<sub>3</sub>, 40% Glycerol

### Product overview

**Description** Anti-BRD9 rabbit polyclonal antibody  
**Applications** ELISA, IHC  
**Immunogen** Fusion protein of human BRD9  
**Reactivity** Human, Mouse  
**Content** 0.7 mg/ml  
**Host species** Rabbit  
**Ig class** Immunogen-specific rabbit IgG  
**Purification** Antigen affinity purification

### Target information

**Symbol** BRD9  
**Full name** bromodomain containing 9  
**Synonyms** PRO9856; LAVS3040  
**Swissprot** Q9H8M2

### Target Background

May play a role in chromatin remodeling and regulation of transcription.

订购热线: 4008-898-798


#### Applications

##### Immunohistochemistry

Predicted cell location: Cell membrane

Positive control: Human thyroid cancer

Recommended dilution: 50-300



The image on the left is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using ml124272(BRD9 Antibody) at dilution 1/80, on the right is treated with fusion protein. (Original magnification:  $\times 200$ )

##### ELISA

Recommended dilution: 5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: [mlbio\\_cn@yeah.net](mailto:mlbio_cn@yeah.net)

网址: [www.mlbio.cn](http://www.mlbio.cn)