

## Anti-CSF1 antibody

<b>Cat. No.</b>	ml263065
<b>Package</b>	25 µl/100 µl/200 µl
<b>Storage</b>	-20°C, pH7.4 PBS, 0.05% NaN <sub>3</sub> , 40% Glycerol

### Product overview

<b>Description</b>	Anti-CSF1 rabbit polyclonal antibody
<b>Applications</b>	ELISA, WB
<b>Immunogen</b>	Synthetic peptide of human CSF1
<b>Reactivity</b>	Human, Rat
<b>Content</b>	0.9 mg/ml
<b>Host species</b>	Rabbit
<b>Ig class</b>	Immunogen-specific rabbit IgG
<b>Purification</b>	Antigen affinity purification

### Target information

<b>Symbol</b>	CSF1
<b>Full name</b>	colony stimulating factor 1 (macrophage)
<b>Synonyms</b>	MCSF; CSF-1
<b>Swissprot</b>	P09603

### Target Background

The protein encoded by this gene is a cytokine that controls the production, differentiation, and function of macrophages. The active form of the protein is found extracellularly as a disulfide-linked homodimer, and is thought to be produced by proteolytic cleavage of membrane-bound precursors. The encoded protein may be involved in development of the placenta. Alternate splicing results in multiple transcript variants.

订购热线: 4008-898-798

**Applications**

**Western blotting**

Predicted band size:60 kDa

Positive control:Human heart tissue lysate

Recommended dilution: 200-1000

Gel: 6%SDS-PAGE


Lysate: 40  $\mu$ g

Lane: Human heart tissue lysate

Primary antibody: ml263065(CSF1 Antibody) at dilution 1/200

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 30 seconds



**ELISA**

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio\_cn@yeah.net

网址: www.mlbio.cn