

Anti-BMP2 antibody

Cat. No.	ml263283
Package	25 µl/100 µl/200 µl
Storage	-20°C, pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol

Product overview

Description	Anti-BMP2 rabbit polyclonal antibody
Applications	ELISA, WB
Immunogen	Synthetic peptide of human BMP2
Reactivity	Human
Content	1.4 mg/ml
Host species	Rabbit
Ig class	Immunogen-specific rabbit IgG
Purification	Antigen affinity purification

Target information

Symbol	BMP2
Full name	bone morphogenetic protein 2
Synonyms	BDA2; BMP2A
Swissprot	P12643

Target Background

This gene encodes a member of the transforming growth factor-beta (TGFB) superfamily. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer, which induces bone and cartilage formation.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:45 kDa

Positive control:Human fetal liver tissue lysate

Recommended dilution: 500-2000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane: Human fetal liver tissue lysate

Primary antibody: ml263283(BMP2 Antibody) at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 20 seconds



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn