

Anti-FAM98A antibody

Cat. No. ml263883
Package 25 µl/100 µl/200 µl
Storage -20°C, pH7.4 PBS, 0.05% NaN₃, 40% Glycerol

Product overview

Description Anti-FAM98A rabbit polyclonal antibody
Applications ELISA, WB
Immunogen Synthetic peptide of human FAM98A
Reactivity Human, Mouse, Rat
Content 0.96 mg/ml
Host species Rabbit
Ig class Immunogen-specific rabbit IgG
Purification Antigen affinity purification

Target information

Symbol FAM98A
Full name family with sequence similarity 98 member A
Synonyms
Swissprot Q8NCA5

Target Background

FAM98A belongs to the FAM98 family. The function of FAM98A remains unknown.

订购热线: 4008-898-798

Applications

Western blotting

Predicted band size:55 kDa

Positive control:Human lung adenocarcinoma tissue, Human liver cancer tissue and Human breast cancer tissue lysates

Recommended dilution: 200-1000

Gel: 8%SDS-PAGE


Lysate: 40 μ g

Lane 1-3: Human lung adenocarcinoma tissue, Human liver cancer tissue and Human breast cancer tissue lysates

Primary antibody: ml263883(FAM98A Antibody) at dilution 1/800

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 10 minutes



ELISA

Recommended dilution:5000-10000

联系电话: 4008-898-798, 021-61725725

联系QQ: 2881505695, 2881505696

邮箱: mlbio_cn@yeah.net

网址: www.mlbio.cn