

兔抗 ACAD8 多克隆抗体

中文名称：兔抗 ACAD8 多克隆抗体

英文名称： Anti-ACAD8 rabbit polyclonal antibody

别名： ARC42; ACAD-8

抗原： ACAD8

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the acyl-CoA dehydrogenase family of enzymes that catalyze the dehydrogenation of acyl-CoA derivatives in the metabolism of fatty acids or branch chained amino acids. The encoded protein is a mitochondrial enzyme that functions in catabolism of the branched-chain amino acid valine. Defects in this gene are the cause of isobutyryl-CoA dehydrogenase deficiency.
Applications:	ELISA, IHC
Name of antibody:	ACAD8
Immunogen:	Fusion protein of human ACAD8

Full name:	acyl-CoA dehydrogenase family, member 8
Synonyms :	ARC42; ACAD-8
SwissProt:	Q9UKU7
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	15-50

