

兔抗 ADAMDEC1 多克隆抗体

中文名称：兔抗 ADAMDEC1 多克隆抗体

英文名称：Anti-ADAMDEC1 rabbit polyclonal antibody

别名：M12.219

相关类别：一抗

抗原：ADAMDEC1

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This encoded protein is thought to be a secreted protein belonging to the disintegrin metalloproteinase family. Its expression is upregulated during dendritic cells maturation. This protein may play an important role in dendritic cell function and their interactions with germinal center T cells.
Applications:	ELISA, IHC
Name of antibody:	ADAMDEC1
Immunogen:	Fusion protein of human ADAMDEC1
Full name:	ADAM-like, decysin 1
Synonyms :	M12.219

SwissProt:	O15204
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and human thyroid cancer
IHC Recommend dilution:	25-100

