

兔抗 ADAMTS6 多克隆抗体

中文名称：兔抗 ADAMTS6 多克隆抗体

英文名称： Anti-ADAMTS6 rabbit polyclonal antibody

别名： ADAM-TS6; ADAMTS-6; ADAM-TS 6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ADAMTS6

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) protein family. Members of the family share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains.
Applications:	ELISA, IHC
Name of antibody:	ADAMTS6
Immunogen:	Synthetic peptide of human ADAMTS6

Full name:	ADAM metalloproteinase with thrombospondin type 1 motif, 6
Synonyms:	ADAM-TS6; ADAMTS-6; ADAM-TS 6
SwissProt:	Q9UKP5
IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	15-50

