

## 兔抗 ADAMTS15 多克隆抗体

中文名称：兔抗 ADAMTS15 多克隆抗体

英文名称： Anti-ADAMTS15 rabbit polyclonal antibody

抗 原： ADAMTS15

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motifs) protein family. ADAMTS family members share several distinct protein modules, including a propeptide region, a metalloproteinase domain, a disintegrin-like domain, and a thrombospondin type 1 (TS) motif. Individual members of this family differ in the number of C-terminal TS motifs, and some have unique C-terminal domains. The protein encoded by this gene has a high sequence similarity to the protein encoded by ADAMTS1 and ADAMTS8.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ADAMTS15

<b>Immunogen:</b>	Fusion protein of human ADAMTS15
<b>Full name:</b>	ADAM metalloproteinase with thrombospondin type 1 motif, 15
<b>SwissProt:</b>	Q8TE58
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human ovarian cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	50-200

