

## 兔抗 ADAMTSL1 多克隆抗体

中文名称：兔抗 ADAMTSL1 多克隆抗体

英文名称：Anti-ADAMTSL1 rabbit polyclonal antibody

别名：C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1

相关类别：一抗

储存：冷冻 (-20℃)

抗原：ADAMTSL1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a secreted protein and member of the ADAMTS (a disintegrin and metalloproteinase with thrombospondin motif) family. This protein lacks the metalloproteinase and disintegrin-like domains, which are typical of the ADAMTS family, but contains other ADAMTS domains, including the thrombospondin type 1 motif. This protein may have important functions in the extracellular matrix. Alternative splicing results in multiple transcript variants encoding distinct proteins.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ADAMTSL1

<b>Immunogen:</b>	Fusion protein of human ADAMTSL1
<b>Full name:</b>	ADAMTS-like 1
<b>Synonyms :</b>	C9orf94; PUNCTIN; ADAMTSR1; ADAMTSL-1
<b>SwissProt:</b>	Q8N6G6
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human lung cancer
<b>IHC Recommend dilution:</b>	15-50

