

兔抗 ADAP1 多克隆抗体

中文名称：兔抗 ADAP1 多克隆抗体

英文名称：Anti-ADAP1 rabbit polyclonal antibody

别名：ArfGAP with dual PH domains 1; GCS1L; CENTA1; p42IP4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ADAP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|--|
| Background: | GTPase-activating protein for the ADP ribosylation factor family (Probable). Binds phosphatidylinositol 3,4,5-trisphosphate (PtdInsP3) and inositol 1,3,4,5-tetrakisphosphate (InsP4). |
| Applications: | ELISA, IHC |
| Name of antibody: | ADAP1 |
| Immunogen: | Fusion protein of human ADAP1 |
| Full name: | ArfGAP with dual PH domains 1 |
| Synonyms: | GCS1L; CENTA1; p42IP4 |
| SwissProt: | O75689 |
| ELISA Recommended dilution: | 5000-10000 |

| | |
|--------------------------------|---|
| IHC positive control: | Human liver cancer and Human thyroid cancer |
| IHC Recommend dilution: | 25-100 |



