

兔抗 ADARB1 多克隆抗体

中文名称：兔抗 ADARB1 多克隆抗体

英文名称： Anti-ADARB1 rabbit polyclonal antibody

别名： RED1; ADAR2; DRABA2; DRADA2

相关类别： 一抗

储存 冷冻（-20℃） 避光

抗原： ADARB1

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型 Unconjugate

技术规格

Background:	This gene encodes the enzyme responsible for pre-mRNA editing of the glutamate receptor subunit B by site-specific deamination of adenosines. Studies in rat found that this enzyme acted on its own pre-mRNA molecules to convert an AA dinucleotide to an AI dinucleotide which resulted in a new splice site. Alternative splicing of this gene results in several transcript variants, some of which have been characterized by the presence or absence of an ALU cassette insert and a short or long C-terminal region.
Applications:	WB
Name of antibody:	ADARB1

Immunogen:	Synthesized peptide derived from internal of human ADAR B1.
Full name:	adenosine deaminase, RNA-specific, B1
Synonyms :	RED1; ADAR2; DRABA2; DRADA2
SwissProt:	P78563
WB Predicted band size:	81 kDa
WB Positive control:	HepG2 cells lysate
WB Recommended dilution:	500-3000

