

兔抗 ADCK1 多克隆抗体

中文名称：兔抗 ADCK1 多克隆抗体

英文名称：Anti-ADCK1 rabbit polyclonal antibody

相关类别：一抗

储存：冷冻（-20℃）

抗原：ADCK1

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).
Applications:	ELISA, IHC
Name of antibody:	ADCK1
Immunogen:	Fusion protein of human ADCK1
Full name:	aaF domain containing kinase 1
SwissProt:	Q86TW2
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and human thyroid cancer
IHC Recommend dilution:	25-100

