

兔抗 ADCK2 多克隆抗体

中文名称：兔抗 ADCK2 多克隆抗体

英文名称： Anti-ADCK2 rabbit polyclonal antibody

别名： AARF

相关类别： 一抗

抗原： ADCK2

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The function of this protein is not yet clear. It is not known if it has protein kinase activity and what type of substrate it would phosphorylate (Ser, Thr or Tyr).
Applications:	WB
Name of antibody:	ADCK2
Immunogen:	Synthesized peptide derived from internal of human ADCK2.
Full name:	aarF domain containing kinase 2
Synonyms :	AARF
SwissProt:	Q7Z695

WB Predicted band size:	69 kDa
WB Positive control:	Jurkat cells and K562 cells lysates
WB Recommended dilution:	500-3000

