

兔抗 ADCY1 多克隆抗体

中文名称：兔抗 ADCY1 多克隆抗体

英文名称：Anti-ADCY1 rabbit polyclonal antibody

别名：AC1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ADCY1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a form of adenylate cyclase expressed in brain. A similar protein in mouse is involved in pattern formation of the brain.
Applications:	ELISA, IHC
Name of antibody:	ADCY1
Immunogen:	Synthetic peptide of human ADCY1
Full name:	adenylate cyclase 1 (brain)
Synonyms :	AC1
SwissProt:	Q08828
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human colon cancer and Human liver cancer

IHC Recommend dilution:

100-300

