

兔抗 ADCY2 多克隆抗体

中文名称：兔抗 ADCY2 多克隆抗体

英文名称： Anti-ADCY2 rabbit polyclonal antibody

别名： AC2; HBAC2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ADCY2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the family of adenylate cyclases, which are membrane-associated enzymes that catalyze the formation of the secondary messenger cyclic adenosine monophosphate (cAMP). This enzyme is insensitive to Ca(2+)/calmodulin, and is stimulated by the G protein beta and gamma subunit complex.
Applications:	ELISA, WB
Name of antibody:	ADCY2
Immunogen:	Synthetic peptide of human ADCY2
Full name:	adenylate cyclase 2 (brain)

Synonyms:	AC2; HBAC2
SwissProt:	Q08462
ELISA Recommended dilution:	5000-10000
WB Predicted band size:	124 kDa
WB Positive control:	293T cell lysate
WB Recommended dilution:	200-1000

