

兔抗 ADCY9 多克隆抗体

中文名称：兔抗 ADCY9 多克隆抗体

英文名称：Anti-ADCY9 rabbit polyclonal antibody

别名：AC9

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ADCY9

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Adenylate cyclase is a membrane bound enzyme that catalyses the formation of cyclic AMP from ATP. It is regulated by a family of G protein-coupled receptors, protein kinases, and calcium. The type 9 adenylyl cyclase is a widely distributed adenylyl cyclase, and it is stimulated by beta-adrenergic receptor activation but is insensitive to forskolin, calcium, and somatostatin.
Applications:	WB, IHC
Name of antibody:	ADCY9
Immunogen:	Synthesized peptide derived from N-terminal of human ADCY9.

Full name:	adenylate cyclase 9
Synonyms :	AC9
SwissProt:	O60503
IHC positive control:	Human lung carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	151 kDa
WB Positive control:	Huvec cells lysate
WB Recommended dilution:	500-3000



