

## 兔抗 ADCYAP1 多克隆抗体

中文名称：兔抗 ADCYAP1 多克隆抗体

英文名称： Anti-ADCYAP1 rabbit polyclonal antibody

别名： adenylate cyclase activating polypeptide 1; PACAP

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ADCYAP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a secreted proprotein that is further processed into multiple mature peptides. These peptides stimulate adenylate cyclase and increase cyclic adenosine monophosphate (cAMP) levels, resulting in the transcriptional activation of target genes. The products of this gene are key mediators of neuroendocrine stress responses. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ADCYAP1
<b>Immunogen:</b>	Synthetic peptide of human ADCYAP1
<b>Full name:</b>	adenylate cyclase activating polypeptide 1

<b>Synonyms:</b>	PACAP
<b>SwissProt:</b>	P18509
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer
<b>IHC Recommend dilution:</b>	20-100
<b>WB Predicted band size:</b>	19 kDa
<b>WB Positive control:</b>	Human cerebrum tissue lysate
<b>WB Recommended dilution:</b>	200-1000



