

兔抗 ADGRB2 多克隆抗体

中文名称：兔抗 ADGRB2 多克隆抗体

英文名称：Anti-BAI2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：ADGRB2

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a a seven-span transmembrane protein that is thought to be a member of the secretin receptor family. The encoded protein is a brain-specific inhibitor of angiogenesis. The mature peptide may be further cleaved into additional products (PMID:20367554). Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	ADGRB2
Immunogen:	Synthetic peptide of human ADGRB2
Full name:	adhesion G protein-coupled receptor B2
Synonyms:	BAI2
SwissProt:	O60241

IHC positive control:	Human gastric cancer and human breast cancer
IHC Recommend dilution:	25-100

