

## 兔抗 ADGRG2 多克隆抗体

中文名称：兔抗 ADGRG2 多克隆抗体

英文名称： Anti-ADGRG2 rabbit polyclonal antibody

别名： adhesion G protein-coupled receptor G2; HE6; EDDM6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ADGRG2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the G protein-coupled receptor family described as an epididymis-specific transmembrane protein. The encoded protein may be proteolytically processed as it contains a motif shown to be a protein scission motif in some members of this family (PMID: 11973329). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011]
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ADGRG2
<b>Immunogen:</b>	Synthetic peptide of human ADGRG2

<b>Full name:</b>	adhesion G protein-coupled receptor G2
<b>Synonyms:</b>	HE6; EDDM6; GPR64; CBAVDX; TM7LN2
<b>SwissProt:</b>	Q8IZP9
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	50-200

