

兔抗 ADGRG6 多克隆抗体

中文名称: 兔抗 ADGRG6 多克隆抗体

英文名称: Anti-ADGRG6 rabbit polyclonal antibody

别名: adhesion G protein-coupled receptor G6; APG1; DREG; VIGR; LCCS9; PR126; GPR126; PS1TP2

储存: 冷冻 (-20℃)

抗原: ADGRG6

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene, which is upregulated in human umbilical vein endothelial cells, encodes a G protein-coupled receptor. Variations in this gene can affect a person's stature. Multiple transcript variants encoding different proteins have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ADGRG6
Immunogen:	Synthetic peptide of human ADGRG6
Full name:	adhesion G protein-coupled receptor G6
Synonyms:	APG1; DREG; VIGR; LCCS9; PR126; GPR126; PS1TP2

SwissProt:	Q86SQ4
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human ovarian cancer
IHC Recommend dilution:	25-100



