

兔抗 ADH4 多克隆抗体

中文名称：兔抗 ADH4 多克隆抗体

英文名称：Anti-ADH4 rabbit polyclonal antibody

别名：ADH-2

相关类别：一抗

抗原：ADH4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes class II alcohol dehydrogenase 4 pi subunit, which is a member of the alcohol dehydrogenase family. Members of this enzyme family metabolize a wide variety of substrates, including ethanol, retinol, other aliphatic alcohols, hydroxysteroids, and lipid peroxidation products. Class II alcohol dehydrogenase is a homodimer composed of 2 pi subunits. It exhibits a high activity for oxidation of long-chain aliphatic alcohols and aromatic alcohols and is less sensitive to pyrazole. This gene is localized to chromosome 4 in the cluster of alcohol dehydrogenase genes.

Applications:	ELISA, IHC
Name of antibody:	ADH4
Immunogen:	Fusion protein of human ADH4
Full name:	alcohol dehydrogenase 4 (class II), pi polypeptide
Synonyms :	ADH-2
SwissProt:	P08319
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human lung cancer
IHC Recommend dilution:	25-100

