

兔抗 ADIPOR1 多克隆抗体

中文名称：兔抗 ADIPOR1 多克隆抗体

英文名称： Anti-ADIPOR1 rabbit polyclonal antibody

别名： CGI45; PAQR1; ACDCR1; CGI-45; TESBP1A

储存： 冷冻（-20℃）

抗原： ADIPOR1

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein which acts as a receptor for adiponectin, a hormone secreted by adipocytes which regulates fatty acid catabolism and glucose levels. Binding of adiponectin to the encoded protein results in activation of an AMP-activated kinase signaling pathway which affects levels of fatty acid oxidation and insulin sensitivity. A pseudogene of this gene is located on chromosome 14. Multiple alternatively spliced transcript variants have been found for this gene.
Applications:	ELISA, IHC
Name of antibody:	ADIPOR1

Immunogen:	Synthetic peptide of human ADIPOR1
Full name:	adiponectin receptor 1
Synonyms :	CGI45; PAQR1; ACDCR1; CGI-45; TESBP1A
SwissProt:	Q96A54
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	20-100



