

兔抗 ADM2 多克隆抗体

中文名称：兔抗 ADM2 多克隆抗体

英文名称：Anti-ADM2 rabbit polyclonal antibody

别名：AM2; dJ579N16.4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ADM2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein which is a member of the calcitonin-related hormones. The encoded protein is involved in maintaining homeostasis in many tissues, acting via CRLR/RAMP receptor (calcitonin receptor-like receptor/receptor activity-modifying protein) complexes. Multiple alternatively spliced variants, encoding the same protein, have been identified.
Applications:	ELISA, IHC
Name of antibody:	ADM2
Immunogen:	Synthetic peptide of human ADM2
Full name:	adrenomedullin 2
Synonyms :	AM2; dJ579N16.4

SwissProt:	Q7Z4H4
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain and Human liver cancer
IHC Recommend dilution:	25-100

