

兔抗 ADORA3 多克隆抗体

中文名称: 兔抗 ADORA3 多克隆抗体

英文名称: Anti-ADORA3 rabbit polyclonal antibody

别名: A3AR

储存: 冷冻 (-20℃)

抗原: ADORA3

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene encodes a protein that belongs to the family of adenosine receptors, which are G-protein-coupled receptors that are involved in a variety of intracellular signaling pathways and physiological functions. The receptor encoded by this gene mediates a sustained cardioprotective function during cardiac ischemia, it is involved in the inhibition of neutrophil degranulation in neutrophil-mediated tissue injury, it has been implicated in both neuroprotective and neurodegenerative effects, and it may also mediate both cell proliferation and cell death. Alternative splicing results in multiple transcript variants. This g

	ene shares its 5' terminal exon with some transcripts from overlapping GeneID:57413, which encodes an immunoglobulin domain-containing protein.
Applications:	ELISA, WB, IHC
Name of antibody:	ADORA3
Immunogen:	Fusion protein of human ADORA3
Full name:	adenosine A3 receptor
Synonyms:	A3AR
SwissProt:	P0DMS8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	36 kDa
WB Positive control:	Mouse liver tissue lysate
WB Recommended dilution:	500-2000

