

## 兔抗 ADRA2B 多克隆抗体

中文名称：兔抗 ADRA2B 多克隆抗体

英文名称： Anti-ADRA2B rabbit polyclonal antibody

别名： ADRA2L1; ADRARL1; ADRA2RL1; ALPHA2BAR; alpha-2BAR

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： ADRA2B

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This intronless gene encodes a seven-pass transmembrane protein. This protein is a member of a subfamily of G protein-coupled receptors that regulate neurotransmitter release from sympathetic nerves and from adrenergic neurons in the central nervous system.
<b>Applications:</b>	IF
<b>Name of antibody:</b>	ADRA2B
<b>Immunogen:</b>	Synthesized peptide derived from internal of human Adrenergic Receptor alpha-2B.
<b>Full name:</b>	adrenoceptor alpha 2B
<b>Synonyms :</b>	ADRA2L1; ADRARL1; ADRA2RL1; ALPHA2BAR; alpha-2BAR

<b>SwissProt:</b>	P18089
<b>IF positive control:</b>	HepG2 cells
<b>IF Recommend dilution:</b>	100-500

