

## 兔抗 ADRB1 多克隆抗体

中文名称：兔抗 ADRB1 多克隆抗体

英文名称：Anti-ADRB1 rabbit polyclonal antibody

别名：RHR, B1AR, ADRB1R, BETA1AR

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ADRB1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The adrenergic receptors (subtypes alpha 1, alpha 2, beta 1, and beta 2) are a prototypic family of guanine nucleotide binding regulatory protein-coupled receptors that mediate the physiological effects of the hormone epinephrine and the neurotransmitter norepinephrine. Specific polymorphisms in this gene have been shown to affect the resting heart rate and can be involved in heart failure.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ADRB1
<b>Immunogen:</b>	Synthetic peptide of human ADRB1
<b>Full name:</b>	Beta-1 adrenergic receptor

<b>Synonyms :</b>	RHR, B1AR, ADRB1R, BETA1AR
<b>SwissProt:</b>	P08588
<b>ELISA Recommended dilution:</b>	1000-5000
<b>IHC positive control:</b>	Human tonsil cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	25-100

