

兔抗 AGAP2 多克隆抗体

中文名称：兔抗 AGAP2 多克隆抗体

英文名称：Anti-AGAP2 rabbit polyclonal antibody

别名：PIKE; GGAP2; CENTG1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：AGAP2

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the centaurin gamma-like family. It mediates anti-apoptotic effects of nerve growth factor by activating nuclear phosphoinositide 3-kinase. It is overexpressed in cancer cells, and promotes cancer cell invasion. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.
Applications:	ELISA, IHC
Name of antibody:	AGAP2
Immunogen:	Synthetic peptide of human AGAP2
Full name:	ArfGAP with GTPase domain, ankyrin repeat and PH d

	omain 2
Synonyms :	PIKE; GGAP2; CENTG1
SwissProt:	Q99490
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human thyroid cancer
IHC Recommend dilution:	50-200

