

兔抗 AGER 多克隆抗体

中文名称：兔抗 AGER 多克隆抗体

英文名称：Anti-AGER rabbit polyclonal antibody

别名：RAGE

相关类别：一抗

储存：冷冻（-20℃）

抗原：AGER

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The advanced glycosylation end product (AGE) receptor encoded by this gene is a member of the immunoglobulin superfamily of cell surface receptors. It is a multiligand receptor, and besides AGE, interacts with other molecules implicated in homeostasis, development, and inflammation, and certain diseases, such as diabetes and Alzheimer's disease. Many alternatively spliced transcript variants encoding different isoforms, as well as non-protein-coding variants, have been described for this gene.
Applications:	ELISA, WB, IHC

Name of antibody:	AGER
Immunogen:	Fusion protein of human AGER
Full name:	Advanced glycosylation end product-specific receptor
Synonyms :	RAGE
SwissProt:	Q15109
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human breast cancer and Human gastric cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	43 kDa
WB Positive control:	231 cells
WB Recommended dilution:	500-2000



