

## 兔抗 AGL 多克隆抗体

中文名称：兔抗 AGL 多克隆抗体

英文名称： Anti-AGL rabbit polyclonal antibody

别名： GDE

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： AGL

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes the glycogen debrancher enzyme which is involved in glycogen degradation. This enzyme has two independent catalytic activities which occur at different sites on the protein: a 4-alpha-glucotransferase activity and a amylo-1,6-glucosidase activity. Mutations in this gene are associated with glycogen storage disease although a wide range of enzymatic and clinical variability occurs which may be due to tissue-specific alternative splicing. Alternatively spliced transcripts encoding different isoforms have been described.
<b>Applications:</b>	ELISA, WB

<b>Name of antibody:</b>	AGL
<b>Immunogen:</b>	Synthetic peptide of human AGL
<b>Full name:</b>	amylo-alpha-1, 6-glucosidase, 4-alpha-glucanotransferase
<b>Synonyms :</b>	GDE
<b>SwissProt:</b>	P35573
<b>ELISA Recommended dilution:</b>	1000-2000
<b>WB Predicted band size:</b>	175 kDa
<b>WB Positive control:</b>	Human fetal liver tissue
<b>WB Recommended dilution:</b>	200-1000

