

兔抗 AGO1 多克隆抗体

中文名称： 兔抗 AGO1 多克隆抗体

英文名称： Anti-AGO1 rabbit polyclonal antibody

别名： Q99; EIF2C; EIF2C1; GERP95

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

克隆类型： rabbit polyclonal

抗原： AGO1

反应种属： Human, Mouse

标记物： Unconjugate

技术规格

Background:	This gene encodes a member of the Argonaute family of proteins which play a role in RNA interference. The encoded protein is highly basic, and contains a PAZ domain and a PIWI domain. It may interact with dicer1 and play a role in short-interfering-RNA-mediated gene silencing. This gene is located on chromosome 1 in a cluster of closely related family members including argonaute 3, and argonaute 4.
Applications:	ELISA, IHC
Name of antibody:	AGO1
Immunogen:	Synthetic peptide of human AGO1
Full name:	argonaute RISC catalytic component 1

Synonyms :	Q99; EIF2C; EIF2C1; GERP95
SwissProt:	Q9UL18
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human cervical cancer and Human brain
IHC Recommend dilution:	50-200

